BORT Helix S Spiraldynamik USG-Bandage

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Med

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Spiraldynamik
intelligent movement

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
REF 053 200

Illustrationen
Figures

1

5

2

6

3

7

4

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	08
FR	Français	Mode d'emploi	12
ES	Español	Instrucciones de uso	16
IT	Italiano	Istruzioni per l'uso	20
NL	Nederlands	Gebruiksaanwijzing	24
CS	Čeština	Návod k použití	28
ET	Eesti	Kasutusjuhend	32
PL	Polski	Instrukcja użytkowania	36
RO	Română	Instrucțiunile de utilizare	40

PDF: ga.bort.com

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BORT Helix S Spiraldynamik Lower Ankle Support

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a support for soft ankle tissue compression consisting of elastic material with additional functional elements. They serve to treat the stabilisation of the upper and lower ankle joint.

Indications

Chronic, post-traumatic or post-operative irritations of soft tissue in the area of lower ankle joint and the sole of the foot, e.g. in case of arthrosis, ligament weakness, foot problems due to static anomalies and defects in the shape of the foot due to dysfunction of the lower ankle joint, e.g. skew pes planotransversus.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks / Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product

- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

- 1 Slip into the dressing aid with your foot.
- 2 Release the support strap. To ease putting the support on, pull it over the dressing aid on the foot. Check the correct fit.
- 3 Then, take hold of the dressing aid at the tip and pull it forwards out of the support. To simplify putting the rest on, we recommend – if possible – placing the foot on the thigh on the opposite side.
- 4 Guide the elastic strap around the ankle joint and over the instep to the outer surface of the foot.
- 5 Guide the strap underneath the sole of the foot to the inner surface of the foot. If necessary, the heel bone can be pulled towards the middle of the body using the other hand in this phase to support the restraints.
- 6 Close the elastic strap by placing the Velcro tip on the velour field above the inner ankle pad. Check and, if necessary, make final corrections to the fit and the individual strap tension of the aid.

Removal

To remove the support, please open the elastic strap, grasp the lower edge of the support and pull it down over the foot.

Tip: Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition

Polyamide (PA), elastodiene / latex (LA), viscose (CV).

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer

Do not iron Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life / Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

MD Medical device |  Single patient – multiple use

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