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The BORT Foot Levator Orthotic Device supports the lifting of the foot (dorsal flexion) in the upper ankle joint when walking. Whilst in the momentum phase, the Foot Levator Orthotic Device supports the resetting of the foot and a physiological movement pattern is facilitated. Both of the partially-elastic bracing elements (A, B) lift the forefoot. These bracing elements are fastened in a figure-of-eight form above the ankle. The figure-of-eight closure additionally stabilises the ankle. The non-elastic bracing element (C) prevents the foot from twisting in the lower ankle area; it lifts the outer edge of the foot.

The permanent elastic, soft and textile structure of the BORT Foot Levator Orthotic Device permits self-perception through the plantar surface of the sole of the foot whilst under strain when walking. The dynamics of walking is possible without any restriction.

Instructions

Applying the BORT Foot Levator Orthotic Device:


The Foot Levator Orthotic Device can be worn directly on the skin or over a thin, closely fitting stocking. It is recommended to apply the Foot levator orthotic device whilst seated.

- Open all Velcro closures on the Foot Levator Orthotic Device and pull it over your foot.
- Reach into the finger loop that is attached to the upper finishing edge on the rear side of the Foot Levator Orthotic Device and pull the Foot Levator Orthotic Device upwards over the heel.

Closing the closures and applying the bracing elements:


Please ensure that the Foot Levator Orthotic Device is correctly positioned on the foot before closing the Velcro bracing elements. The Foot Levator Orthotic Device must be positioned correctly on the heel. For an efficient running of the bracing elements, the seams along the back of the foot as well as in the heel area must run centrally.

Important: When applying the Foot Levator Orthotic Device, the heel of the foot is to remain on the floor for the entire application procedure.

Step 1

Place the foot on the floor and now place bracing elements A, B and C onto the floor as shown in the drawing.


Step 2

Now close the Velcro closures located on the lower leg in the following order:

- First close the lowermost Velcro point.
- Then close the uppermost Velcro point.
- Finally close the Velcro point in the centre.

Adjust the tension of the Velcro closures of the Foot Levator Orthotic Device on the lower section if required in order to achieve a comfortable yet secure hold of the Foot Levator Orthotic Device.

Step 3

Push the thigh or knee as far forward as possible. An angle from the foot to the lower leg amounting to between 20 and 30° is ideal.

First guide bracing element A upwards in a vertical manner, then over the back of the foot, above the medial malleolus and around the lower leg whilst applying tension. Then secure the bracing element onto the hook element attached to the rear of the brace.

Step 4

Fasten bracing element B in the opposite direction to bracing element A. First guide bracing element B upwards in a vertical manner, then over the back of the foot, above the lateral malleolus and around the lower leg whilst applying tension. Then fasten bracing element B onto bracing element A via the hook element attached to the rear of the brace. Secure bracing element B on the hook element of hook A with as much alignment as possible in order to prevent stockings from becoming damaged.

Step 5

Finally, guide the non-elastic bracing element C around the outer side of the lower leg whilst applying upward tension. In doing so, slightly lift the outer edge of the foot with the bracing element. Then fasten bracing element C to the velour Velcro area of the Foot Levator Orthotic Device. **Important:** bracing element C must run upwards in a vertical manner. It stabilises the ankle and should prevent the ankle from twisting outwards whilst walking.

Indications: Peroneal paralysis of any severity, especially for mild to medium paralysis of the tibialis anterior muscle and the flexor hallucis longus (degree of severity 2 – 4).

Contraindications:

Peripheral arterial occlusive disease (PAOD), lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid, sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body, latex allergy.

Composition: Knitted fabric: 65 % polyamide, 28 % elastodies, 7 % viscose

Important Note: No reutilisation – this medical device is meant for the care of one patient only. In case of doctor’s order to wear support during night, avoid hindering blood circulation. In case of numbness loosen support and remove it if necessary. If complaints continue, please contact the doctor. While wearing the support do not use creams or ointments. This can damage the material.

The product contains latex and can trigger allergic reactions.

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