# BORT KubiFX light Ellenbogen-Orthese

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**REF** 055 350

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D05535012023-051006 ML I Rev. 01

### DE BORT KubiFX light Ellenbogen-Orthese

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### EN BORT KubiFX Light elbow brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a brace for immobilisation of the elbow joint in a bent position.

#### Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of proximal sulcus ulnaris syndrome, descending tendinopathy, following bursectomy.

#### Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks / Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- do not close it too tightly during longer periods of rest, e.g. when sleeping
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

#### Fitting instructions

Before putting the BORT KubiFX Light Elbow Brace on for the first time, both internal aluminium splints must be formed in accordance with the desired joint position. Following this, check the positioning of the elbow pad and guide the fixing straps through the guide loops up to the anti-slip lock. The brace can then be attached to the arm. The tip of the elbow (olecranon) must lie on the integrated bend (curve) of the aluminium splints. The brace's fixing straps are then closed with the desired tension. Both middle belt straps in the vicinity of the elbow bends cross distally to the elbow bend in the lower arm region, thus not causing any undesired pressure on the nerves and blood vessels located there. Finally, it must be checked that, the elbow joint is fixed securely, but no blood flow congestion is incurred.

#### Removal

To remove the product, loosen the belt straps and take the brace off the arm.

#### Material composition

PUR foam (PUR), cotton (CO), polyamide (PA)  
The sewn in textile label provides the precise material composition.

#### Cleaning information

Delicate wash 30 °C Do not bleach (only wash with colour or mild detergents) Do not dry in a tumble dryer Do not iron  
 Do not clean chemically

Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life / Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 05.2023

Medical device | Single patient – multiple use

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