BORT KubiFX lang Ellenbogen-Orthese

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REF 055 400

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D055400I2022-051004 ML I Rev. 01

DE BORT KubiFX lang Ellenbogen-Orthese

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EN BORT KubiFX Long elbow brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for immobilisation of the elbow joint with immobilisation of the proximal radial ulnar joint.

Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of proximal sulcus ulnaris syndrome, descending tendinopathy, epicondylitis, following bursectomy.

Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks / Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- do not close it too tightly during longer periods of rest, e.g. when sleeping
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient

Fitting instructions

The brace is single-piece and has a freely positionable elbow pad which is attached to the brace using Velcro. Before putting the BORT KubiFX Elbow Brace on, the two pre-bent aluminium splints located in the splint bags on the reverse side of the brace, must be formed in accordance with the desired bending position in which the elbow joint is to be fixed. The pre-moulded aluminium splint integrated into the brace on the lower (palm) side of the arm and can also be reshaped if needed. Following this, the elbow pad is positioned and the fixing straps, which are equipped with a lock against slipping through, are guided through the D rings.

The brace can then be attached to the arm, the tip of the elbow must lie on the bend (curve) of the aluminium splints on the pad. The first Velcro fastener is guided over the commissure between the thumb and the index finger and attached to the brace using Velcro. The D rings already guided using fixing straps are tightened and closes, starting on the forearm.

Finally, check that the elbow joint is fixed securely, but no blood flow congestion is incurred. The reversible strap attached on the inner and outer

lining in the vicinity of the antecubital fossa serves as an additional stretching inhibitor and is adjusted individually. In order to ease attachment, the upper arm slip-on aid can be removed by cutting it out if required.


Removal

To remove the product, loosen the belt straps and take the brace off the arm.

Material composition

PUR foam (PUR), cotton (CO), polyamide (PA)
The sewn in textile label provides the precise material composition.

Cleaning information

 Delicate wash 30 °C  Do not bleach (only wash with colour or mild detergents)  Do not dry in a tumble dryer  Do not iron  Do not clean chemically
Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life / Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 05.2022

Medical device |  Single patient – multiple use

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