# BORT Air Walker Achillo

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**Med**  
Mitten im  
gesunden  
Leben

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
**REF 100 320**

Illustrationen  
Figures


1


4

2


5

3

6

## Sprachen Languages

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	04
<b>EN</b>	English	Instructions for use .....	07
<b>FR</b>	Français	Mode d'emploi .....	10
<b>ES</b>	Español	Instrucciones de uso .....	13
<b>IT</b>	Italiano	Istruzioni per l'uso .....	16
<b>NL</b>	Nederlands	Gebruiksaanwijzing .....	19
<b>CS</b>	Čeština	Návod k použití .....	22
<b>ET</b>	Eesti	Kasutusjuhend .....	25
<b>PL</b>	Polski	Instrukcja użytkowania .....	28
<b>RO</b>	Română	Instrucțiuni de utilizare .....	31

 PDF: [ga.bort.com](https://ga.bort.com)

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

**Intended purpose**

This medical device is a lower leg foot brace for complete immobilisation of the lower leg and foot region in a defined, adjustable position.

**Indications**

Achilles tendon rupture, pre-operative or post-operative/post-traumatic immobilisation of the lower leg and foot region, post-operative in case of Achilles tendon joint.

**Contraindications**

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

**Application risks/Important notes **

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- the Walker is appropriate for a maximum body weight of 130 kg
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

**Fitting instructions**

Before going through steps 1–6, pull one of the two pairs of understockings provided over your foot/lower leg.

**Information:**

The opaque, see-through protective caps are only used for transporting the Walker. They must be removed prior to use. When ready, the protective caps no longer need be attached to the blue pumps at the side. They could well prove to be an obstacle when walking.

- 1 Open the circumferential Velcro fasteners and remove the front cover.
- 2 Position the wedges provided to set the therapeutically necessary equinus position in the footrest. One wedge segment is equivalent to an elevation of approx. 10 mm. Not all wedge segments must be used. For final positioning, the individual segments can be fixed using the attached adhesive points. The wedge height can be adapted to the course of therapy.
- 3 Open the black padding liner and place the foot and lower leg in the Walker. Pay attention to good heel positioning. The black padding liner is then closed, pay attention it sitting crease-free.
- 4 Now, adapt the grey front cover by attaching it centrally to the black padding liner again using Velcro. Close the three circumferential Velcro fasteners. Start with the fastener on the ankle joint, then the fastener over the instep area and finally the lower leg fastener. If necessary, re-tighten again.
- 5 Now, fill the integrated air chambers using the blue pumps at the side. The air chambers can be filled individually. Start with the (medial) pump, then the outer (lateral) pump.
- 6 If you wish to reduce the air pressure, you can release air. For this, fold the grey protective cap next to the pumps upwards and press the silver-coloured button.

**Accessories available:**

Ref. 100 370 Understocking

Ref. 100 380 Protective cover against wetness

Ref. 100 390 Protective cover for night-time

**Removal**

If you wish to remove the Walker, fold the grey protective cap next to the pumps upwards and press the silver-coloured button. Then, open the circumferential Velcro fasteners and remove the front cover. Finally, open the black padding liner and remove your foot from the Walker.

**Tip:** Close the Velcro fastener in order to avoid early wear to the Velcro surface.

**Material composition**

Polyester (PES), polyurethane (PUR), polyamide (PA)

The sewn in textile label provides the precise material composition.

**Cleaning information**

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Hand wash the shoe using lukewarm water and mild detergents.

**Guarantee**

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

**Useful life/Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

**Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

**Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

**Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 07.2023

**MD** Medical device |  Single patient – multiple use

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