# BORT MalleoStabil Air/Gel-Orthese

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
**REF** 100 540

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D100540I2023-10I007 ML I Rev. 01

### DE BORT MalleoStabil Air/Gel-Orthese

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### EN BORT MalleoStabil Air/Gel Brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is an ankle brace for the stabilisation of the upper ankle joint on one plane. Its function is to restrict pronation and supination via two connected stabilisation elements. The stabilisation elements are attached using adjustable belts.

#### Indications

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament ruptures, ankle fractures (Weber A), chronic ligament instability, shortening of time spent in plaster, ligament injuries, upper ankle joint distortion, fibular ligament rupture, pronation and/or supination trauma.

#### Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

- 1 When putting the brace on, ensure that the lateral shells are precisely centred with the lower leg. Please set the strap for the sole of the foot to fit the patient's circumstances in order to achieve effective stabilisation. Lay the wide lower strap circumferentially around the shells and attach it with Velcro.
- 2 Lay the upper strap circumferentially around the shells and attach it with Velcro. If required, adjust the lower strap so that the brace is firmly in position.
- 3 When the swelling is reduced, it is possible that the MalleoStabil Air/Gel Brace becomes loose. If required, adjust it.

- 4 Use of the air/gel component  
Remove the gel component from the bottom upwards. When re-attaching the gel component, observe that it is placed directly at its original location. Its correct placing ensures that the plastic edges are covered.

#### For cold therapy

- Remove the air/gel component and put it into a deep-freezer for an hour.
- **Information:** The gel doesn't harden.
- Re-attach the air/gel component at its original location in the brace.
- Do not bring the air/gel component into direct contact with the skin.
- Always wear a sock for best possible comfort.

#### Removal

You only need to open the two circumferential belts to remove the BORT MalleoStabil Brace.

**Tipp:** Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

#### Material composition

80 % Polypropylene (PP), 20 % polyamide velour (PA)

#### Cleaning information

Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  
 Do not iron  Do not clean chemically  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2023

Medical device | Single patient – multiple use

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