BORT MalleoXpress Active

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
REF 100 580

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D10058012023-101006 ML | Rev. 01

DE BORT MalleoXpress Active

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EN BORT MalleoXpress Active

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is an ankle joint brace for the immobilisation of the ankle joint in a defined position.

Indications

Post-operative protection after ligament suture, conservative treatment of ankle joint distortions and ligament ruptures, immobilisation in the case of inflammatory irritations, arthritis, ligament injuries, supination/pronation trauma.

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks / Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

- 1 Open the circumferential Velcro fasteners and release the Boa fastener by pulling it outwards until it „clicks“ audibly. The fastener is now unlocked and the pulley can be drawn out. Additionally, release the metatarsal belt and open the wide fastener over the metatarsus.

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- 2 Place your foot in the fully opened support from the front. Firstly, close the wide fastening over the metatarsus. Pay attention to the support sitting crease-free and the ankle joint being placed at an angle of 90°.

- 3 Following this, the two circumferential straps are closed on the upper edge of the brace. Firstly, the rear strap, then the front strap.

- 4 Following this, guide the metatarsal belt inwards over the metatarsus and hook it in there.

- 5 Push the Boa fastener wheel towards the body with an audible „click“. The fastener is now engaged and the brace can be closed. Now turn the Boa fastener wheel clockwise in order to achieve stabilisation of your ankle joint. The more you tighten the Boa system, the better the stabilisation. Please pay attention that, if you fasten it too tightly, neither pressure sores, erythema nor influence of the blood circulation are caused.

Removal

- 6 To remove the brace, the Boa fastener must be released. For this purpose, release the Boa fastener by pulling it outwards until it „clicks“ audibly. You can then draw the pulley out of the Boa casing with the textile strap and loosen the MalleoXpress Active brace.

- 7/8 Finally open the remaining fastenings, release the metatarsal belt and remove the brace from the leg.

Information: The plastic stabilisation elements positioned at the side of the brace are not thermoplastically remouldable/mouldable.

Material composition

Polyamide (PA), elastane (EL)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use

of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

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Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2023

 Medical device |  Single patient – multiple use

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