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EN

Immobilisation is a basic principle in the treatment of injuries and overstrain since it greatly accelerates the recovery process and prevents secondary complications. Hard-frame knee braces support damaged tissue structures (joints, ligaments, capsules, cartilage), stabilise, guide and hold, thereby contributing to a positive recovery process.

Instruction for use:

Before fitting the brace please set the extension and flexion locks on the medial and lateral tooth segment joint. The extension is preset to 0°.

A different extension or flexion stop can be set by loosening the pre-inserted mounting screw, removing the limitation wedge and securing with the mounting screw.

Then open all the fixing straps and fit the brace with the leg slightly flexed. Position the brace such that the rotation point of the brace joints is in the horizontal plane at the height of the middle of the patella and in the vertical plane along the Nietert compromise axis (ratio 40 : 60).

First adjust the frontal fastener on the thigh such that the brace cannot slide backwards. Then close the frontal infrapatellar fastener on the lower leg. Close the fixing straps – one by one - from bottom to top.

Adjustment of metal frame and tibia clasp:

The brace can be adapted to the patient’s physiognomy as desired. This is done by bending the metal frame and heating the tibia clasp which is thermoplastic at 130 °C with a hot air gun. *When shaping/bending the aluminium frame, it is essential to ensure that the mechanical joints are parallel in all planes, otherwise the function of the mechanical joints may be impaired.*

Indications:

Serious and / or complex instability of the knee joint, meniscus injuries, genu recurvatum, functional pre-and / or post-operative treatment of ligaments ruptures (ACL, PCL, collateral ligaments).

Contraindications:

Peripheral arterial occlusive disease (PAOD), lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid, sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body.

Composition:

60 % aluminium, 22 % nylon, 18 % polyurethane

Important Note:

No reutilization – this medical device is meant for the care of one patient only.

In case of doctor’s order to wear brace during night, avoid hindering blood circulation. In case of numbness loosen brace and remove if necessary. If complaints continue, please contact the doctor. If pressure points, abrasions or similar occur, have your orthopaedic technician or doctor check the fit of the brace immediately.

Safety information:

Under the terms of the German Medical Devices Act the manufacturer’s product liability lapses when the product is re-used and/or processed by the service provider.

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