BORT GenuXpress

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Med
Mitten im
gesunden
Leben

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

REF 100 700

Illustrationen
Figures


1


2


3


4

5

6

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	08
FR	Français	Mode d'emploi	12
ES	Español	Instrucciones de uso	16
IT	Italiano	Istruzioni per l'uso	20
NL	Nederlands	Gebruiksaanwijzing	24
CS	Čeština	Návod k použití	28
ET	Eesti	Kasutusjuhend	32
PL	Polski	Instrukcja użytkowania	36
RO	Română	Instrucțiuni de utilizare	40

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace in a frame construction for guiding and stabilising the knee joint using a 4-point principle with belt straps and adjustable joints for limiting the scope of movement.

Indications

Serious and/or complex instability of the knee joint, meniscus injuries, genu recurvatum, functional pre-and/or post-operative treatment of ligaments ruptures (ACL, PCL, collateral ligaments).

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

Fitting instructions

Instructions for use for patients:

The brace can be worn directly on the skin or over the clothes. Pay attention to avoid creasing the clothes in the vicinity of the brace as this can result in pressure sores.

- 1** Open all the brace's fasteners. You can either do this by opening the quick-release fasteners or by unlooping the harnesses.
- 2** Sit down and slip forwards to the edge of the chair. Bend your knee joint slightly (approx. 30°).

3 Place the brace on your leg. Ensure correct positioning and even surface contact. The joints should be positioned at the height of the centre of the patella (centre of the kneecap).

4 Close the fasteners in the following order, always from the centre of the knee downwards, then from the centre of the knee upwards:

1. Front fastener close to the knee on the lower leg
2. Front fastener close to the knee on the thigh
3. Rear fastener close to the knee on the lower leg
4. Fastener away from the knee on the lower leg
5. Fastener close to the knee on the thigh
6. Fastener away from the knee on the thigh

Whilst tightening the fastener, smoothen it in the tension direction with your free hand. This avoids twisting the brace on the leg.

Information: Should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician. Allow yourself time to become accustomed to the brace and increase the wearing time gradually. If you have questions regarding the use and wearing time, please refer to your physician or technician.

Removal

To remove the aid, open the harnesses and then take the brace off the leg.

Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

Information for technicians:

Important note: The brace may only be adapted by expert staff.

Adaption of the BORT GenuXpress knee brace on the patient:

Adaption of the BORT GenuXpress knee brace on the patient:

1 Open all the brace's fasteners. You can either do this by opening the quick-release fasteners or by unlooping the harness. The strap length can be set by shortening the straps.

2 The patient should sit down and slip forwards to the edge of the chair, the knee joint is bent approx. 30°.

3 Place the brace on the patient's leg. Ensure correct positioning according to Nietert. The joints should be positioned at the height of the centre of the patella. If required, the side condyle pads (10 mm material thickness) can be replaced with a thinner version (7 mm material thickness) in order to regulate the clear span in the knee joint area.

4 Close the fasteners in the following order, always from the centre of the knee downwards, then from the centre of the knee upwards:

1. Front fastener close to the knee on the lower leg
2. Front fastener close to the knee on the thigh
3. Rear fastener close to the knee on the lower leg
4. Fastener away from the knee on the lower leg
5. Fastener close to the knee on the thigh
6. Fastener away from the knee on the thigh

Whilst tightening the fastener, smoothen it in the tension direction with your free hand. This avoids twisting the brace on the leg.

Important: Let the patient try walking. Check for possible pressure sores, if required the brace frame can be reshaped. Be careful not to damage the brace frame

5+6 Setting the extension/flexion limit:
If required, it can be necessary to restrict the scope of the brace's movement.

1. Remove the outer joint cover. For this, place both thumb tips on the upper edge of the joint cover **and push the joint cover away downwards whilst applying slight pressure.**
2. The joints are marked with "A" and "B". Remove the existing limit stop from the joint and take the required limit stop from the limit stop sets provided.
They are also marked with the letters "A" and "B", the angular degrees noted on the limit stop wedges. Extension wedges are marked with an "E", flexion wedges with an "F".
3. Insert the limit stop into the joint.
4. Finally, re-attach the joint cover, place the cover on the joint for this purpose. Press the cover on evenly so that the release flap protrudes slightly upwards and slide the cover upwards until it engages.

Available spare part: REF 991120 Quick-release fastener

Important note

Please always insert an extension and flexion limit stop into the joints. For unrestricted or free movement, insert the 0° extension and 140° flexion limit stop into the joints.

Use of the brace without inserted limit stops results in damages to the joint and dysfunction of the brace. Ensure the use of identical extension and flexion stops in joint A and B in order to avoid an uneven load of the joint structure.

The following restrictions are available:

- Extension: 0°, 5°, 10°, 15°, 20° and 25°
- Flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

Material composition

Frame pad: 60% EVA foam, 40% polyester

Fastening straps: 100% polyamide

Fastening straps pad: 74% polyurethane, 26% polyether

Cleaning information

Can be cleaned using mild detergents and lukewarm water. Close the Velcro fasteners in order to avoid damaging items of clothing, for example. Air dry it and avoid heat.

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron
 Do not clean chemically

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2023

 Medical device |  Single patient – multiple use

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