# BORT ROM Knieorthese

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**Med**  
Mitten im  
gesunden  
Leben

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

**REF 100 900**

Illustrationen  
Figures


[1]


[6]


[2]


[7]


[3]


[8]


[4]


[9a]


[9b]

[5]

[10]

## Sprachen Languages

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	04
<b>EN</b>	English	Instructions for use .....	09
<b>FR</b>	Français	Mode d'emploi .....	14
<b>ES</b>	Español	Instrucciones de uso .....	19
<b>IT</b>	Italiano	Istruzioni per l'uso .....	24
<b>NL</b>	Nederlands	Gebruiksaanwijzing .....	29
<b>CS</b>	Čeština	Návod k použití .....	34
<b>ET</b>	Eesti	Kasutusjuhend .....	38
<b>PL</b>	Polski	Instrukcja użytkowania .....	42
<b>RO</b>	Română	Instrucțiuni de utilizare .....	47

 PDF: [ga.bort.com](https://ga.bort.com)

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### Intended purpose

Brace for the mobilisation of the knee joint.

### Indications

Postoperative, posttraumatic early functional mobilisation with restricted movement of the knee joint, condition after damage to bone/capsule/ligament structures of the knee joint.

### Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

### Instructions for use for patients

The brace can be worn directly on the skin or over the clothes. Avoid easing the clothes in the vicinity of the brace, this may result in pressure sores.

#### Step 1

[1], [2] Open all the brace's fasteners. This can either be done by opening the blue quick action fasteners or by unlooping the harnesses.

#### Step 2

[3] Sit down and slip forwards to the edge of the chair. Bend your knee joint slightly.

**Step 3**

**[4]** Place your leg in the pre-set and adapted BORT ROM knee brace, adaptations/re-adjustments must be made by expert staff.

Ensure correct positioning and even surface contact. The joints should be positioned on the side of the knee joint at the height of the centre of the patella (middle of the kneecap).

**Step 4**

**[4]** Firstly, close the wide fastener beneath the knee joint in order to fix the brace on the leg.

**Step 5**

**[5]** Close the fasteners in the following order, always starting from the centre of the knee:

1. Fastener close to the knee on the lower leg (#1)
2. Fastener close to the knee on the thigh (#2)
3. Fastener away from the knee on the lower leg (#3)
4. Fastener away from the knee on the thigh (#4)

Whilst tightening the fastener, smoothen it in the tension direction with your free hand. This avoids twisting the brace on the leg.

**Information:** should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician

Allow yourself time to become accustomed to the brace and increase the wearing time gradually.

If you have questions regarding the use and wearing time, please refer to your physician or technician.

**Removal**

To remove the aid, open the harnesses and then take the brace off the leg. Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

**Information for the technician**

Important note: The brace may only be adapted by expert staff.

**Adaption of the BORT ROM knee brace on the patient****Step 1**

**[1], [2]** Open all the brace's fasteners. This can either be done by opening the blue quick action fasteners or by unlooping the harnesses. The strap length can be set shortening the straps. For this, the crocodile fasteners can be removed at the respective end of the strap. After successful setting of the fastener lengths, they can be fixed using Velcro tape.

**Step 2**

**[3]** The patient should sit down and slip forwards to the edge of the chair, the knee joint is bent slightly.

### Step 3

[6], [7] Initially set the required brace length. For this, first open the lock on the medial and lateral joint splints on the thigh and lower leg. The length can be set on the side splints in the “Unlock” position. The splint should end 5-7 cm before the crotch medially on the thigh. **Important:** Secure the length set by locking in the “Lock” position.

### Step 4

[4] Ensure correct positioning and even surface contact. The joints should be positioned on the side of the knee joint at the height of the centre of the patella (middle of the kneecap).

In order to maintain a form fit of the joint splints on the leg, the side joint splints can be bent up to the contour of the leg.

### Step 5

#### [8] Use as a rest splint:

For positioning, the brace can also be fully blocked in a defined position. For this, slide the blue slider directly on the joint to the “Unlock” position, set the required extent of bending and finally lock the joint again. Immobilisation can be made at  $-10^{\circ}$ ,  $0^{\circ}$ ,  $10^{\circ}$ ,  $20^{\circ}$  or  $30^{\circ}$ .

#### [10] Use as a ROM splint:

The treating physician provides information regarding any necessary restrictions. Slide the blue slider directly on the joint to the “Unlock” position to free the joints completely. Then, follow step 6 to set the scope of movement

### Step 6

[9a], [9b] The therapeutically necessary extension/flexion restriction is made directly on the joint.

For this, the blue bar on the limit stops is folded upwards.

[10] Simply push the limit stop downwards to unlock it and move it to the desired angular degree. The possible degree settings are printed directly on the joint. Please observe that the blue lever must again be folded downwards to lock the angle adjustment.

### Step 7

[4] Firstly, close the wide fastener beneath the knee joint in order to fix the brace on the leg.

### Step 8

[5] Close the fasteners in the following order, always starting from the centre of the knee:

1. Fastener close to the knee on the lower leg (#1)
2. Fastener close to the knee on the thigh (#2)
3. Fastener away from the knee on the lower leg (#3)
4. Fastener away from the knee on the thigh (#4)

Whilst tightening the fastener, smoothen it in the tension direction with your free hand. This avoids twisting the brace on the leg.

**Important:** If possible, let the patient try walking. Check for possible pressure sores, if required the brace frame can be reshaped. Be careful not to damage the brace frame.

### **Step 9**

Then, instruct the patient in the use of the aid.

### **Material composition**

Frame pad: 60% EVA foam, 40% polyester

Fastening straps: 100% polyamide

### **Cleaning information**

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

### **Guarantee**

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

### **Useful life/Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

**Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 02.2024

**MD** Medical device |  Single patient – multiple use

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