# BORT X-Walker kurz

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**Med**  
Mitten im  
gesunden  
Leben

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
**REF 101 350**

Illustrationen  
Figures


1


2


3


4

5

6

## Sprachen Languages

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	04
<b>EN</b>	English	Instructions for use .....	08
<b>FR</b>	Français	Mode d'emploi .....	12
<b>ES</b>	Español	Instrucciones de uso .....	16
<b>IT</b>	Italiano	Istruzioni per l'uso .....	20
<b>NL</b>	Nederlands	Gebruiksaanwijzing .....	24
<b>CS</b>	Čeština	Návod k použití .....	28
<b>ET</b>	Eesti	Kasutusjuhend .....	32
<b>PL</b>	Polski	Instrukcja użytkowania .....	36
<b>RO</b>	Română	Instrucțiuni de utilizare .....	40

 PDF: [ga.bort.com](http://ga.bort.com)

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a lower leg foot brace for complete immobilisation of the lower leg and foot region in a specified position.

### **Indications**

Pre-operative or post-operative/post-traumatic immobilisation of the forefoot and metatarsus region, metatarsal fracture, soft tissue injuries, distortions.

### **Contraindications**

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

### **Application risks/Important notes** ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- the Walker is appropriate for a maximum body weight of 130 kg
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

### **Fitting instructions**

**1** Open the circumferential Velcro fasteners and remove the front cover. In order to make it easier to put on later, we recommend twisting the Velcro straps downwards, thus easing application of the brace.

**2** Open the padding liner and place the foot and lower leg in the Walker. Pay attention to good heel positioning, it should be placed securely in the very back of the heel rest. The padding liner is

then closed on the foot and in the lower leg region using the Velcro fastener, pay attention to it sitting crease-free.

### Information for the technician:

#### Setting the foot length:

If required, the foot length can be adapted to the patient's circumstances by moving the front binding. For this purpose, the front binding is raised and pulled upwards gently to disengage it. Following this, it can be moved and be locked in the desired position by pressing it downwards gently.

**3** Then, adapt the grey front cover again centrally on the padding liner. The valve at the upper end should be freely accessible.

Close the four circumferential Velcro fasteners in the order specified below:

1. Start with the fastener on the ankle joint
2. Close the lower leg fastener
3. Close the strap in the metatarsus/toe region.

If necessary, re-tighten again.

**4+5** If necessary, remove the protective cap from the valve and slide the blue valve slider downwards. Place the vacuum pump provided on the valve. In order to draw the air out of the liner and generate the vacuum, now pump several times until the pump no longer jumps back. Finally, remove the pump from the valve and re-attach the protective cap.

You should now feel surface stabilisation in the lower leg and foot region. If necessary, re-tighten the fastening straps again.

**Tip:** To ensure a good fit, we recommend that you pay attention to even distribution of the filling, which consists of a large number of small balls. By fluffing up the padding liner before putting the brace on or smoothing out the filling with your hand, you can achieve even distribution of the filling and optimise its fit.

#### Accessories available:

REF 100 370 Understocking

REF 100 380 Protective cover against wetness

REF 100 390 Protective cover for night-time

#### Removal

**6** To remove the aid, remove the protective cap from the valve. To reduce the vacuum, pull the blue valve slider upwards. The air then flows into the vacuum chamber audibly.

Then, open the harnesses, remove the front cover and open the liner to remove the brace from the leg.

Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

**Information:** The vacuum cushion can be removed from the liner cover for cleaning purposes. To remove the liner cover, remove the liner from the X-Walker. Open the zip on the reverse side completely. Guide the valve through the opening in the cover and remove the vacuum cushion from the liner cover. Before cleaning, close the zip again.

After cleaning and thorough drying of the liner cover, open the zip and insert the vacuum cushion into the cover. Guide the valve through the opening in the liner cover.

Attention: Check the valve position before inserting it. The valve must point in the direction of the valve passage in the liner cover.

### Material composition

Brace rest:	100% polypropylene (PP)
Vacuum chamber:	100% polyvinylchloride (PVC)
Vacuum chamber filling:	100% polystyrene (EPS)
Liner (cover):	15% polyester (PES), 55% cotton (CO), 30% PUR foam (PUR)
Fasteners:	100% polyamide (PA)

### Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron  
 Do not clean chemically

Clean the brace rest manually with lukewarm water, wash the liner cover with mild detergents at a maximum of 30 °C. Close the Velcro fasteners in order to avoid damaging items of clothing, for example. Air dry whilst avoiding the impact of heat.

Before cleaning the liner cover, ensure removal of the vacuum cushion.

### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the responsible national body in your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact)

**Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

**Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 06.2024

**MD** Medical device |  Single patient – multiple use

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