BORT OA-Xpress

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





Med
Mitten im
gesunden
Leben

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REF 101 600

Illustrationen
Figures

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	08
FR	Français	Mode d'emploi	12
ES	Español	Instrucciones de uso	16
IT	Italiano	Istruzioni per l'uso	20
NL	Nederlands	Gebruiksaanwijzing	24
CS	Čeština	Návod k použití	28
ET	Eesti	Kasutusjuhend	32
PL	Polski	Instrukcja użytkowania	36
RO	Română	Instrucțiuni de utilizare	40

 PDF: ga.bort.com

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Stand: 09.2023

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for relieving and guiding the knee joint via the 3-point principle with belt straps and adjustable joints for restricting the scope of movement.

Indications

Partial relief of the medial or lateral compartment in case of level II-III gonarthrosis, bone bruise (post-traumatic bone marrow edema), medial/lateral gonarthrosis as a result of a tibial plateau fracture.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

Fitting instructions

Instructions for use for patients:

The brace can be worn directly on the skin or over the clothes. Pay attention to avoid creasing the clothes in the vicinity of the brace as this can result in pressure sores.

- 1** Open all the brace's fasteners. You can either do this by opening the quick action fasteners or by unlooping the belt straps.
- 2** Sit down and slip forwards to the edge of the chair. Bend your knee joint slightly (approx. 30°).

3 Place the brace on your leg. Ensure correct positioning. The joint should be positioned at the height of the centre of the patella (middle of the kneecap) and on the side of the leg.

4 Close the fasteners in the following order, always from the centre of the knee downwards, then from the middle of the knee upwards:

1. Fastener close to the knee on the lower leg
2. Fastener away from the knee on the lower leg
3. Fastener close to the knee on the thigh
4. Fastener away from the knee on the thigh

Whilst tightening the fastener, smoothen it in the tension direction with your free hand. This avoids twisting the brace on the leg.

Information: should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician. Allow yourself time to become accustomed to the brace and increase the wearing time gradually. You can remove the brace during longer periods of rest. If you have questions regarding the use and wearing time, please refer to your physician or technician.

5 Tip: You can find our video, which shows the attachment process, on our homepage.

Removal

To remove the aid, open the harnesses and then take the brace off the leg.

Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

Information for the technician:

Important note: the brace may only be adapted by expert staff! The BORT OA-Xpress can provide for axial differences of up to 5°.

Application of model A:

- Use in case of varus dislocation/relief of the medial compartment on the right leg (the brace is positioned laterally).
- Use in case of varus dislocation/relief of the medial compartment on the left leg (the brace is positioned medially).

Application of model B:

- Use in case of varus dislocation/relief of the lateral compartment on the right leg (the brace is positioned medially).
- Use in case of varus dislocation/relief of the medial compartment on the left leg (the brace is positioned laterally).

Adaption of the BORT OA-Xpress knee brace on the patient:

1 Open all the brace's fasteners. You can either do this by opening the quick action fasteners or by unlooping the belt straps. The strap length can be set using the shorter straps.

2 The patient should sit down and slip forwards to the edge of the chair, the knee joint is bent approx. 30°.

3 Place the brace on the patient's leg and use the enclosed key to adapt the brace to the patient's leg axis. Ensure correct positioning. The polycentric joint should be positioned at the height of the centre of the patella and on the side of the leg. If required, the side condyle pad (10 mm material thickness) can be replaced with a thinner version (7 mm material thickness).

4 Close the fasteners in the following order:

1. Proximal fastener on the lower leg
2. Distal fastener on the lower leg
3. Distal fastener on the thigh
4. Proximal fastener on the thigh

Whilst tightening the fastener, smoothen it in the tension direction with your free hand. This avoids twisting the brace on the leg.

5 Tip: You can find our video, which shows the attachment process, on our homepage.

6 Use the enclosed key to set the required correction pressure when the brace is positioned and the patient is sitting. For this purpose, adjustment screws are provided on the rear beneath the joint. Avoid overcorrection!

Important: let the patient try walking. If required, the fitting correction pressure can be readjusted.

7 Setting the extension/flexion limit:

if required, it can be necessary to restrict the scope of the brace's movement.

1. For this, remove the condyle pad
2. Remove the screw for the limit stop to be set
3. Remove the limit stop from the joint and insert the required angular degree
4. Fix the stop used with the screw provided

The following restrictions are available:

- Extension: 0°, 5°, 10°, 15°, 20° and 25°
- Flexion: 0° (180°), 30° 45°, 60°, 75° and 90°

Available spare part: REF 991120 Quick-release fastener

Material composition

Frame pad: 60% EVA foam, 40% polyester

Fastening straps: 100% polyamide

Fastening straps pad: 74% PU foam, 26% polyether

Cleaning information

Can be cleaned using mild detergents and lukewarm water. Close the Velcro fasteners in order to avoid damaging items of clothing, for example. Air dry it and avoid heat.

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 09.2023

 Medical device |  Single patient – multiple use

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