# BORT Dorsale Intrinsic Plus Schiene (D.I.P.S.)

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**REF** 102 500

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D102500|2024-10|007 ML I Rev. 01

### DE BORT Dorsale Intrinsic Plus Schiene (D.I.P.S.)

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### EN BORT Dorsal Intrinsic Plus Splint (D.I.P.S.)

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a wrist brace for finger fixation and immobilisation of the wrist and fingers.

#### Indications

Pre-operative, post-operative, post-traumatic, e.g. metacarpus and proximal finger phalanx in the area of the long fingers (e.g. a boxer fracture), ligament lesions in the long fingers in the proximity of the metacarpophalangeal joints, injuries and illnesses which necessitate immobilisation of the long fingers (e.g. infections).

#### Contraindications

Circulatory disorders with various origins, lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated, tendon injuries, CRPS (Morbus Sudeck), burns, compartment syndrome.

#### Application risks/Important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
  - This brace may only be adapted by trained expert staff. If required, have the brace reshaped or remodelled in accordance with the specialist doctor's instructions. For this, open the zip and remove the aluminium splint from the pad cover.
  - do not close it too tightly during longer periods of rest, e.g. when sleeping
  - in case of numbness, loosen the medical device and remove it if necessary
  - in case of persistent complaints, consult the physician or specialist retailer
  - use the medical device in accordance with therapeutic needs
  - only use other products simultaneously after consultation with your physician
  - do not make any changes to the product
  - do not wear it on open wounds
  - do not use in case of intolerance of one of the materials used
  - no re-use. This medical aid is intended for treating one patient only
  - whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material

#### Fitting instructions

Before putting the brace on, open all the Velcro straps. Place the brace on the back of the hand and close the Velcro fasteners on the wrist and lower arm. Then, close the Velcro fastener which fixes the metacarpus to the dorsal brace. Finally, close the circular Velcro fastener which fixes the fingers to the

dorsal brace. If required, the velour tape provided can additionally be attached in the wrist area with a Velcro fastener.




#### Removal

Open all the Velcro fasteners and remove the brace. Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

#### Material composition

Polyurethane foam (PUR), polyamide (PA)  
The sewn in textile label provides the precise material composition.

#### Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron  Do not clean chemically  
Remove the aluminium splint before washing. Close the Velcro fasteners in order to avoid damaging other items of clothing. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

 Medical device |  Single patient – multiple use

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