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**EN**

Affections of the pelvis ring (e.g. slackness of the connective tissue of several pelvic bones) and in connection with it instability are a frequently occurring after-effect of the reorientation process of the hormonal household of people passing the age of 50. Females - not at east caused by pregnancy - are certainly more frequently affected than males.

### Instructions for use:

Position the brace in the middle on the pelvis so that the upper edge of the brace is not higher than the upper iliac crest or protruding into the waist. Insert a finger in the finger loop provided on the fastener of the abdominal plate. Now pull the fastener evenly forwards.

Press the left half of the fastener onto the abdomen and place the right half of the fastener over the left half until you can attach it on the hook and loop pad on the left fastener half. When securing, withdraw the finger first from the left and then from the right finger loop.

Now grasp both tension straps in both hands and pull them simultaneously and steadily forwards until the desired force effect (stabilisation) is produced. Then secure both tension straps with the hook and loop fastener on the abdominal plate.

**Indications:** Symphysis insufficiency, symphysis rupture, structural loosening in the area of the iliosacral joints/pubic symphysis, after birth.

**Contraindications:** Lymphatic drainage disorders including unclear soft tissue swelling, sensory disturbances and circulatory disorders in the specific body region, skin diseases in the specific part of the body.

**Composition:** 44 % polyamide, 25 % polyester, 25 % viscose, 4 % spandex, 2 % cotton

### Important Note:

No re-use – this medical aid is intended for treating one patient. If the support is to be worn at night, observe that the circulation is not impeded. In case of feelings of numbness, loosen the support and remove it if necessary. In case of persistent complaints, you should consult a doctor.

Neither use any creams nor ointments when wearing the support as they could damage the material.

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