# BORT Narbenbruch-Bandage

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**REF** 104 050

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D10405012022-101005 ML | Rev.01

### DE BORT Narbenbruch-Bandage

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### EN BORT Incisional hernia support

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is an elastic, circumferential brace which is suitable for treating damages caused by a hernial ring in the torso region.

#### Indications

Incisional hernia, incisional insufficiency, abdominal muscle infirmity/paralysis.

#### Contraindications

Unclear soft tissue swellings, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material
- In order to avoid damages to the knitted support fabric, only attach the pad after successfully determining the size and trial fitting. For this, remove the protective label and attach the pad in the support. Goods damaged by the Velcro fastening cannot be returned.

#### Fitting instructions

Attach the pad to the support with Velcro and place it on the hernial ring. The pad helps to reposition (move back) the hernia. For this purpose, the pad is reinforced with an inserted plastic core and can be positioned freely on the entire support thanks to the hook strap. For attaching the aid, the patient can be placed on an appropriate examination table in psaos position: For this, lie down flat on your back and put your legs up. The lower legs are placed at right angles to the thighs on a surface which should be as high as the length of the thighs.

The incisional hernia should then no longer emerge from the abdomen. Place the support around the body and close it.

#### Removal

To remove, open the fastener and then remove the support.

#### Material composition

Polyamide (PA), polyester (PES), elastane (EL)

The sewn in textile label provides the precise material composition.

#### Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Close the Velcro fastening to avoid damage to other laundry items.  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2022

Medical device | Single patient – multiple use

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