BORT Stoma-Bandage

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
REF 104 090

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D10409012023-101006 ML I Rev.01

DE BORT Stoma-Bandage

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EN BORT Stoma Support

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a support, a circumferential compression torso bandage with reinforcement elements for supporting the abdominal wall.

Indications

Paralysis and hernia of the abdominal wall, post-operatively, post ileostomy and colostomy status with concomitant connective tissue weakness

Contraindications

Unclear soft tissue swellings, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks / Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

- 1 If used with a pad, firstly attach it with Velcro on the inner side of the support in the region of the stoma recess.
- 2 Put the support on loosely. Loosen both Velcro fastenings for the stoma recess.
- 3 Position the opening around the stoma.
- 4 Guide the stoma bag through the opening and close both Velcro fastenings again.
- 5 Finally, close the support so that a sufficient, pleasant stabilisation effect is achieved.

6 Double Velcro fastening for easy placement and positioning of the stoma recess, respectively on the left and right in the abdominal region. Available accessory: REF 104 099 Stoma protector

Removal

To remove the BORT Stoma support, open all latching halves in the abdominal region, also open all Velcro fasteners in the vicinity of the stoma covering, then remove the support.

Material composition

Polyamide (PA), viscose (CV), cotton (BW), elastodiene/latex (LA), elastane (ED)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 60 °C Do not bleach Do not dry in a tumble dryer
 Do not iron Do not clean chemically
Remove the reinforcing plate and fix the Velcro closure to avoid damage to other laundry.
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2023

Medical device | Single patient – multiple use

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