
BORT StabiloPlus Rückenbandage

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Med
Mitten im
gesunden
Leben

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REF 104 730

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	03
EN	English	Instructions for use	06
FR	Français	Mode d'emploi	09
ES	Español	Instrucciones de uso	12
IT	Italiano	Istruzioni per l'uso	15
NL	Nederlands	Gebruiksaanwijzing	18
CS	Čeština	Návod k použití	21
ET	Eesti	Kasutusjuhend	24
PL	Polski	Instrukcja użytkowania	27
RO	Română	Instrucțiuni de utilizare	30

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for partial relief due to lordosis reduction of the lumbar spine.

Indications

Moderate lumboschialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade I spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history), i.e. after intervertebral disc surgery with mild to moderate prolapse, spinal canal stenosis, instability syndrome.

Contraindications

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body, latex allergy.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

When delivered, the BORT StabiloPlus Back support is pre-set for maximum stabilisation. The bridging segment is attached to the support via the cloth bag, the more rigid version of the splints can be found in the splint bags. A splint set which is more flexible is included in the packaging. Open the support, loosen the side tension straps, place the support around the waist and attach both closing sides using Velcro evenly one above the other with a slight pull of the hand straps. Ensure that the support is placed centrally and the lower edge of the support ends in the inguinal region (forward bend). This is the only means to achieve an optimum effect from the BORT StabiloPlus Back support. Tighten the two lower belt straps evenly first; as a next step tighten the upper belt straps, enabling a guaranteed comfortable fit of the BORT StabiloPlus. If overtightened, the support may become wrinkled. During the course of therapy, it may

become necessary to reduce the stabilising effect of the BORT StabiloPlus Back support. This decision is taken by the treating physician. For this, remove the more rigid version of the splints from the externally attached splint bag by folding the restraints back and removing the splints through the slit in the splint bag. The shorter, more flexible splints can then be inserted. The bridging segment can also be removed from its cloth bag and replaced with two splints. For this, remove the cloth bag including the bridging segment from the BORT StabiloPlus Back support, take the bridging segment out of the cloth bag from the rear slit and insert two of the longer, more flexible splints into the holders in the cloth bag designated for this purpose. Finally, attach the cloth bag back to the BORT StabiloPlus; here the wide side facing the direction of the head.

Removal

To remove the brace, open the tension straps and the front fastening and take the brace off.

Material composition

Elastodiene/latex (LA), polyamide (PA), viscose (CV)

The sewn in textile label provides the precise material composition.

 The product contains latex and can trigger allergic reactions.

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron
 Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items. Before cleaning, remove the reinforcement element from the brace.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

MD Medical device |  Single patient – multiple use

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