# BORT Stablio ISG-Orthese

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
**REF** 104 750

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D104750|2023-07|006 ML I Rev. 01

1

### DE BORT Stablio ISG-Orthese

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### EN BORT Stablio sacroiliac brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a function-securing brace for stabilising the pelvis and the sacroiliac joint.

#### Indications

Sacroiliac joint arthrosis, sacroiliac joint instability, sacroiliac joint blockage, myalgias and tendopathies in the pelvic region, structural abnormalities after lumbar spinal spondylodoses.

#### Contraindications

Unclear soft tissue swellings, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

#### Application risks/Important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
  - remove it during longer periods of rest, e.g. when sleeping
  - in case of numbness, loosen the medical device and remove it if necessary
  - in case of persistent complaints, consult the physician or specialist retailer
  - use the medical device in accordance with therapeutic needs
  - only use other products simultaneously after consultation with your physician
  - do not make any changes to the product
  - do not wear it on open wounds
  - do not use in case of intolerance of one of the materials used
  - no re-use. This medical aid is intended for treating one patient only
  - whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

Open the brace and loosen the tension straps. If required, you can use the Velcro fastener to adapt the position of the pads to the patient's individual circumstances.

1 The wide side of the pads points upwards here, the narrow side in the direction of the posterior.

Then place the opened brace circumferentially around the hip. Regarding the correct height, position the aid centrally on the pelvis so that the upper edge (sewn-in label) of the brace is not higher than the upper iliac crest or extends to the waist.

Push your fingers into the finger flaps intended for this purpose at the ventral plate fastening. Then use both hands to pull the latch plates forwards evenly. Place the left latching half on the stomach and position the right latching half far enough over the left one so that you can Velcro it onto the left latching half. When attaching it using Velcro, remove the fingers first from the left and then the right finger flap.

Then take both tension straps in both hands and pull them forwards

simultaneously and evenly until the desired physical influence (stabilisation) has been built-up. Following this, attach both tension straps firmly to the ventral plate using Velcro.

#### Removal

To remove the brace, open the tension straps and the front fastening and take the brace off.

#### Material composition

Polyamide (PA), cotton (BW), elastodiene/latex (LA)  
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

#### Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items. Remove pads prior to the washing procedure, if required clean the brace with a damp cloth.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).  
You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 07.2023

Medical device | Single patient – multiple use

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