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The rear bands support the damaged tissue structures. The compression provided by the orthosis aids the supporting function of the ligaments, thus relieving the tension on the spine. The internal stabilising structure, acting in combination with the circular band straightens lordosis, thus simultaneously avoiding hyperlordosis of the lumbar spine. It immobilises the painful section of the spine, relieves muscular tension and aids the healing process. Thanks to the “modular principle” the stabilising effect of the orthosis can be adapted during the course of treatment. Degeneration of muscles and loss of strength can be prevented thanks to the activating effect of this elastic compression orthosis.

Instruction for use:

At the time of delivery the BORT StabiloPlus stabilising orthosis is set to maximum stabilisation. The stabilising frame is secured to the orthosis via the fabric pockets; the special pockets hold the more rigid version of the splints. A set of more flexible splints is also provided in the package. Position the orthosis around the waist and join the two fasteners together with the Velcro, apply slight traction with the aid of the guide rings. Ensure that the orthosis is positioned centrally on the torso and that the lower edge of the orthosis, in the groin area (the folds of the thigh), is closed. This is the only way to obtain the full effect of the BORT StabiloPlus stabilising orthosis. First of all fasten the two lower straps then tighten the upper straps to guarantee a comfortable fit. Fastening too tightly may cause creases in the orthosis. To unfasten the orthosis first loosen all the straps then open fully. During treatment it may be necessary to reduce the tightness of the BORT StabiloPlus stabilising orthosis.

To do this, remove the more rigid version of the splints from the outside pockets by pulling the straps back and taking out the splints through the special openings in the pocket. Now the shorter, more flexible splints can be inserted.

The transition element and the abdominal reinforcing plate can be removed, if necessary, by loosening the band's Velcro fasteners, thus enabling the orthopedic specialist to shape it thermoplastically.

Indications:

Severe lower back pain with muscular weakness and significant hernia protrusion/prolapse, without absolute indication for surgery, severe lumbar radicular or pseudo-radicular syndrome (refractory to conservative treatment), grade 2 spondylolisthesis with recurring lower back pain, severe lumbar deformity with facet syndrome/arthrosis, post-surgery for herniated disks without resulting paralysis with removal of the hernia, decompression of the spinal canal without deformity, single level herniated disk surgery, short stenosis of the spinal canal without significant paralysis, vertebral fractures resulting from minor trauma, superior end plate deformities in the lumbar spine.

Contraindications:

Lymphatic drainage disorders including unclear soft tissue swelling, sensory disturbances and circulatory disorders in the specific body region, skin diseases in the specific part of the body, latex allergy.

Composition:

57 % polyamide, 39 % elastodien, 4 % viscose

Important Note:

No re-use – this medical aid is intended for treating one patient. If the support is to be worn at night, observe that the circulation is not impeded. In case of feelings of numbness, loosen the support and remove it if necessary. In case of persistent complaints, you should consult a doctor. Neither use any creams nor ointments when wearing the support as they could damage the material.

The product contains latex and can trigger allergic reactions.

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