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BORT StabiloPlus Bridging Brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for relief due to lordosis reduction and functional mobilisation of the lumbar spine.

Indications

Severe lumboischialgia with muscular failure and substantial disc protrusion/prolapse with no clear OP indication, severe radicular, pseudo-radicular lumbar syndrome (resistant to conservative therapy), grade II spondylolisthesis with recurrent lumboischialgia, severe lumbar deformity with facet syndrome/arthrosis, state after thoracic surgery with no remaining paresis in case of cleared out thorax, single-level, without deformity, single-level thoracic surgery, spinal canal stenosis with a short walking distance without noteworthy paresis, vertebral fractures of a limited extent, superior end plate deformities in the lumbar spine.

Contraindications

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body, latex allergy.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

When delivered, the BORT StabiloPlus bridging brace is pre-set for maximum stabilisation. The bridging frame is attached to the support via the cloth bag, the more rigid version of the splints can be found in the splint bags. A splint set which is more flexible is included in the packaging. The additional abdominal pressure plate is also attached to the support. If required, it can be removed. It serves to increase abdominal counterpressure.

Open the support, loosen the side tension straps, place the support around the waist and attach both closing sides using Velcro evenly one above the other with a slight pull of the hand straps. Ensure that the support is placed centrally on the torso and that the lower support edge (antecubital fold on the thigh) ends in the inguinal region. This is the only means to achieve an optimum effect from the BORT StabiloPlus bridging brace. Tighten the two lower belt straps evenly first; as a next step tighten the upper belt straps, enabling a guaranteed comfortable fit of the brace. If overtightened, the support may become wrinkled. During the course of therapy, it may become necessary to reduce the stabilising effect of the BORT StabiloPlus bridging brace. This decision is taken by the treating physician.

For this, remove the more rigid version of the splints from the externally attached splint bag by folding the restraints back and removing the splints through the slit in the splint bag. The shorter, more flexible splints can then be inserted.

The bridging segment and the abdominal reinforcement plate can also be removed from the support, as required, by loosening the Velcro fastenings and be individually thermoplastically adapted by the orthopaedic technician.

Removal

To remove the brace, open the tension straps and the front fastening and take the brace off.

Information for the technician:

If required, the option exists to reduce the overall size of the abdominal ventral plate by cutting it to size via the circumferential recess (groove). The cutting edge should then be cleaned. Additionally, the abdominal ventral plate can be shaped thermally at the edge by heating it to approx. 160°C.

Material composition

Polyamide (PA), elastodiene/latex (LA), viscose (CV)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron

Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items. Before cleaning, remove the reinforcement element from the brace.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for

the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2023

MD Medical device |  Single patient – multiple use

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