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EN BORT Wrist Support with Hemmed Thumb Opening

(LATEX) The product contains latex and can trigger allergic reactions.

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose
This medical device is a wrist support for compression of soft tissue.

Indications
Minor distortions, strains, after removal of plaster cast, tenosynovitis.

Contraindications
Circulatory disorders with various origins, lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes
The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
– remove the medical device prior to radiological examinations
– remove it during longer periods of rest, e.g. when sleeping
– in case of numbness, loosen the medical device and remove it if necessary
– in case of persistent complaints, consult the physician or specialist retailer
– use the medical device in accordance with therapeutic needs
– only use other products simultaneously after consultation with your physician
– do not make any changes to the product
– do not wear it on open wounds
– do not use in case of intolerance of one of the materials used
– no re-use – this medical aid is intended for treating one patient
– whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions
Place the support on the metacarpophalangeal joint affected and guide the thumb through the intended thumb opening at the same time.

Removal
To remove the support, grasp the fabric in the palm area and pull the support off over the fingers.

Material composition
Polyamide (PA), elastodiene/latex (LA), viscose (CV)
The sewn in textile label provides the precise material composition.

Cleaning information
 Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer
 Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

(MD) Medical device | Single patient – multiple use

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