# BORT SellaFix N Mittelhand-Daumenschiene

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**REF** 112 040

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D112040I2024-10I006 ML I Rev.01

### DE BORT SellaFix N Mittelhand-Daumenschiene

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### EN BORT SellaFix N Metacarpus Thumb Splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb.

#### Indications

Post-traumatic, post-operative conditions, i.e. after sprain, contusion, ligament lesions (skier's thumb), immobilisation in the case of rheumatic-inflammatory processes (rhizarthrosis).

#### Contraindications

Circulatory disorders with various origins, lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- **Important:** Please do not subject the SellaFix N metacarpus thumb splint to direct heat or temperatures exceeding 50 °C (e.g. hot water, vehicles, the sauna, heating) as this may result in damaging the material
- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions for the technician/physician:

Please observe that the correction options provided for the SellaFix® thumb brace should be made by expert staff or a physician. To optimise the fit, the brace can either be formed to fit the patient's circumstances via targeted heating to 100 °C for approx. 2 minutes, e.g. with a hot air dryer (please observe sufficient distance of approx. 15 cm here) or using full surface heating with an ordinary commercial heating oven. The brace may not be heated in a microwave oven!

Open the brace before heating it. After being heated, the brace, which is flexible in the meantime, can be directly formed on the skin after cooling down briefly. After fitting the heated brace, close the fasteners, starting at the metacarpophalangeal joint distally towards the thumb and model the plastic shell to the patient's anatomical circumstances. The Velcro fasteners guarantee stepless adjustability.

#### Fitting instructions for the user:

Slip the brace over the wrist and position the thumb in the thumb tube. Now close the circumferential Velcro fastener on the wrist and the thumb. The Velcro fastener guarantees stepless adjustability.

#### Removal

To remove the brace, open the circumferential Velcro on the thumb tube and wrist and slip the brace off in the direction of the fingers.

#### Material composition

59% polyamide (PA), 26% polyurethane (PUR), 15% PE foam

#### Cleaning instruction

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer  
 Do not iron Do not clean chemically  
Do not use fabric conditioner. Dry in the fresh air.  
Close the Velcro fastening to avoid damage to other laundry items.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-kontakt](http://www.bort.com/md-eu-kontakt).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

Medical device | Single patient – multiple use

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