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### **EN** BORT Wrist Support with Thumb Opening

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### **Intended purpose**

This medical device is a wrist support for compression of soft tissue.

#### **Indications**

Minor distortions, strains, after removal of plaster cast, tenosynovitis.

#### **Contraindications**

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

#### **Application risks/Important notes**

- The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
  - remove it during longer periods of rest, e.g. when sleeping
  - in case of numbness, loosen the medical device and remove it if necessary
  - in case of persistent complaints, consult the physician or specialist retailer
  - use the medical device in accordance with therapeutic needs
  - only use other products simultaneously after consultation with your physician
  - do not make any changes to the product
  - do not wear it on open wounds
  - do not use in case of intolerance of one of the materials used
  - no re-use. This medical aid is intended for treating one patient only
  - whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

#### **Fitting instructions**

Place the support on the metacarpophalangeal joint affected and guide the thumb through the intended thumb opening at the same time. The Stabulo fleece strap can be attached individually dosed. For this, attach the Stabulo fleece strap to the velour on the support intended for this purpose using the narrow Velcro strip. The Stabulo fleece strap can then be twisted around the wrist circumferentially.

#### **Removal**

To remove the support, please loosen the Stabulo fleece strap. Grasp the knitted fabric in the palm area and pull the support off over the fingers. Close the Velcro fastener in order to avoid early wear to the Velcro surface.

#### **Material composition**

Polyamide (PA), elastodiene/latex (LA), viscose (CV)  
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

#### **Cleaning information**

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items.

#### **Guarantee**

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### **Useful life/Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2023

Medical device | Single patient – multiple use

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