# BORT SellaTex® Plus - Quervain

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
**REF** 112 160

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D112160|2024-10|007 ML | Rev. 01

### DE BORT SellaTex® Plus – Quervain

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### EN BORT SellaTex® Plus – Quervain

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a wrist brace with a thumb fixation for immobilisation.

#### Indications

**Tendovaginitis stenosa de Quervain**, pre-operative, post-traumatic, post-operative, e.g. after thumb ligament reconstruction, immobilisation in case of rheumatic-inflammatory processes, in case of degenerative diseases such as rhizarthrosis and carpal tunnel syndrome

#### Contraindications

Circulatory disorders with various origins, lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
  - do not close it too tightly during longer periods of rest, e.g. when sleeping
  - in case of numbness, loosen the medical device and remove it if necessary
  - in case of persistent complaints, consult the physician or specialist retailer
  - use the medical device in accordance with therapeutic needs
  - only use other products simultaneously after consultation with your physician
  - do not make any changes to the product
  - do not wear it on open wounds
  - do not use in case of intolerance of one of the materials used
  - no re-use. This medical aid is intended for treating one patient only
  - whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material

#### Fitting instructions

Open all Velcro fasteners.

1 Slip your hand into the hand brace whilst observing that the thumb is guided through the thumb element. If required, the expert staff can adapt the integrated aluminium splints on the back of the hand, in the palm region and on the thumb element individually. Start by closing the strap (1) around the wrist, followed by the metacarpus strap (2) and the strap on the lower arm (3). Finally, close both straps (4) in the thumb region. Pay attention to firm and comfortable fitting.

2 If downgrading the brace is required for therapeutic reasons, the attached thumb element can be removed. In order to do this, loosen the Velcro flaps on the attached thumb element.

3 Slip the thumb element off over the tip of the thumb.

#### Removal

Open the Velcro fasteners, loosen all the straps and slip the brace off towards the fingers. Then, close the fasteners again in order to avoid early wear to the Velcro surfaces.

#### Material composition

polyamide (PA), Polyurethane (PUR), cotton (BW)  
The sewn in textile label provides the precise material composition.

#### Cleaning information

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Remove the aluminium splints before washing. Close the Velcro fasteners in order to avoid damaging other items of clothing. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

Medical device | Single patient – multiple use

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