# BORT SellaXpress kurz

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**REF** 112 840

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D112840|2024-10|006 ML | Rev.01

### DE BORT SellaXpress kurz

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### EN BORT SellaXpress short

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb.

#### Indications

Post-traumatic, post-operative conditions; e.g. following ligament lesions, sprain; immobilisation in the case of rheumatic-inflammatory processes such as rhizarthrosis.

#### Contraindications

Circulatory disorders with various origins, lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
  - do not close it too tightly during longer periods of rest, e.g. when sleeping
  - in case of numbness, loosen the medical device and remove it if necessary
  - in case of persistent complaints, consult the physician or specialist retailer
  - use the medical device in accordance with therapeutic needs
  - only use other products simultaneously after consultation with your physician
  - do not make any changes to the product
  - do not wear it on open wounds
  - do not use in case of intolerance of one of the materials used
  - no re-use. This medical aid is intended for treating one patient only
  - whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

When delivered, the brace is anatomically pre-formed. To optimise the fitting shape, the integrated aluminium stabilisation element should be bent to the patient's situation by bending it. This should be done exclusively by specialist orthopaedic stores or the physician. For fitting, slip the brace over the wrist and position the thumb in the thumb tube. Ensure that the circumferentially grasping positioning on the thumb contacts evenly without causing pressure sores. Now close the circumferential fastener on the thumb, then the wide fastener on the wrist. The Velcro fastener enables stepless setting of the tightness. As the aluminium stabilisation element is formable, this brace can repeatedly be refitted during the course of therapy in order to retain an optimum fit.

The brace can be used in water without any problems. If it comes into contact with saline water, rinse it afterwards with ample fresh water.

#### Removal

To remove the support, please open the circumferential straps and slip the thumb brace off.

#### Material composition

Polyamide (PA), polyurethane (PUR)

The sewn in textile label provides the precise material composition.

#### Cleaning information

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Do not use fabric conditioner. Dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

Medical device | Single patient – multiple use

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