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EN BORT Osgood-Schlatter Knee Joint Support

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a textile support for facilitating movement of the patella. The pressure pad around the kneecap and the adjustable tension element can achieve correction and relief of the kneecap.

Indications

Influence of the movement of the patella, Osgood-Schlatter disease, attachment tendinosis of the tibial tuberosity ligament, gonarthrosis, arthritis, articular effusion, irritations.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Slip the support over the forward section of the foot and then pull it up to the knee. Here, the upper arc of the silicone half ring should be placed directly above the kneecap, the lower silicone pad then surrounds the shin roughness beneath the knee joint (tibial tuberosity ligament). Following this, attach the correction restraint by gradually moving it upwards to the designated Velcro area.

Removal

To remove the support, loosen the correction restraint, grasp the lower edge of the support on the splints and pull it downwards over the foot. Close

the Velcro fastener in order to avoid damaging other items, e.g. to further clothing.

Material composition

Polyamide (PA), elastodiene/latex (LA), elastane (EL), viscose (CV)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Close the Velcro fastening before washing to avoid damage to other laundry items.

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer
 Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 05.2023

MD Medical device | Single patient – multiple use

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