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#### EN BORT StabiloGen Eco with a Grip Top

#### BORT StabiloGen Eco with a Grip Top short

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a knee support for compression of soft tissues with a local pressure pad around the kneecap consisting of flexible material.

#### Indications

Chronic post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e.g. in case of gonarthrosis, arthritis, ligament instability, meniscopathy.

#### Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid position, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

#### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

1 Fold the grip top outwards, grasp the upper edge of the support.

2 Then pull it over your foot up to the knee with the pad facing forwards.

3 The pad should be placed comfortably over the kneecap. Then, flap the grip top back towards the leg so that the silicone nodules are directly on the skin.

#### Removal

4 To remove the support, please fold the grip top outwards, grasp the lower edge of the support at the side on the splints and pull it downwards over the foot.

#### Material composition

Polyamide (PA), elastodiene/latex (LA), elastane (EL), viscose (CV)  
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

#### Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer

Do not iron Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2023

Medical device | Single patient – multiple use

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