BORT Greiflippe aus Aluminium

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REF 117 170

DE BORT Greiflippe aus Aluminium

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Stand: 07.2023

FR BORT Préhension en aluminium

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This aid is a gripper or a "helping hand" which extends the arm, enabling its user to grip objects.

Indications

Paraplegie, restricted movement.

Contraindications

Muscular failure/Lower arm region or hand muscle dysfunction.

Application risks / Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will instruct you regarding its use.

- do not make any changes to the product
- do not grip in case of open palm wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- be careful when gripping fragile or sensitive items

Operating instructions

Grasp the blue hand-held device and operate the moveable lever using your fingers. If you press it, the gripper closes and if you apply no pressure it opens again. Depending on the surface quality/form, a maximum of some 1.5 kg can be lifted. The BORT aluminium gripper weighs approx. 260 g and can provide an extension of approx. 75 cm.

Removal

Please store the gripper at room temperature.

Status: 07.2023

FR BORT Préhension en aluminium

4-Kant-Rohr, Federn: Aluminium 50%
Grippers: PVC (polyvinyl chloride) 5%

Grip: ABS plastic (acrylonitrile butadiene styrene) 20%
Springs: Iron/zinc 25%

Cleaning information

To clean, wipe with a damp cloth. Do not use any chemical cleaning agents.

Garantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

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