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Vierhoekige stang, veren: aluminium	50%
Grijphulp: polyvinylchloride (PVC)	5%
Handgriep: ABS-kunststof (acrylnitril-butadieen-styrol)	20%
Veren: ijzer / zink	25%

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Indikationen
Paraplegie, restricted movement.

Contraindications
Muscular failure/Lower arm region or hand muscle dysfunction.

Application risks / Important notes
This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will instruct you regarding its use.

- remove the medical device prior to radiological examinations
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not grip in case of open palm wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- be careful when gripping fragile or sensitive items

Useful life/Lifetime of the product
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use.

Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 11.2019

Stand: 11.2019

EN BORT
Aluminium Grippers

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose
This aid is a gripper or a "helping hand" which extends the arm, enabling its user to grip objects.

Indications
Paraplegia, restricted movement.

Contraindications
Muscular failure/Lower arm region or hand muscle dysfunction.

Application risks / Important notes
This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will instruct you regarding its use.

- remove the medical device prior to radiological examinations
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not grip in case of open palm wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- be careful when gripping fragile or sensitive items

Operating instructions
Grasp the blue hand-held device and operate the moveable lever using your fingers. If you press it, the gripper closes and if you apply no pressure it opens again. Depending on the surface quality/form, a maximum of some 1.5 kg can be lifted. The BORT aluminium gripper weighs approx. 260 g and can provide an extension of approx. 75 cm.

Removal
Please store the gripper at room temperature.

Status: 11.2019

Stand: 11.2019

Material composition		
Square pipes, springs: Aluminium	50%	
Grippers: PVC (polyvinyl chloride)	5%	
Grip: ABS plastic (acrylonitrile butadiene styrene)	20%	
Springs: Iron/ zinc	25%	

Cleaning information
To clean, wipe with a damp cloth. Do not use any chemical cleaning agents.

Guarantee
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

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