BORT OmoTex Traction

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REF 120 800

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D12080012023-071007 ML 1 Rev.01

DE BORT OmoTex Traction

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EN BORT OmoTex Traction

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a shoulder joint compression support with additional functional elements. A weight band serves as an additional functional element for extending the subacromial space.

Indications

Impingement syndrome, chronic post-traumatic or post-operative soft tissue irritations in the shoulder joint region, periarthropathia humeroscapularis (PHS)

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- ensure that you remove the medical device prior to radiological examinations, it contains magnetic material
- the weight band incorporates four 250 g weight elements for individual traction weighting, in accordance with the course of therapy these weight elements can be reduced
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

Slip the support over the upper arm and the afflicted shoulder until the humerus head is placed in the recess intended for this purpose. Guide the elastic harness past the front and back of the head and guide the healthy arm through the loop so that the strap runs diagonally over the chest and does not cut into the neck and is positioned beneath the opposite shoulder on the chest. Alternatively, the harness can also be guided with both parts behind the upper body to the opposite shoulder. Slip the health arm through the strap so that, like a rucksack holder, it runs around the opposite shoulder and over the back. Regulate the tension of the elastic harness to the desired extent using the Velcro fasteners. The interior, two-part silicone pad can be positioned freely using the Velcro

fastener. It should support and guide the humerus head respectively from the front and back.

To remove the weight band from the shoulder support, pull the sealing tape out of the loop and remove the Velcro fastener.

Removal

To remove the OmoTex Traction support, open all the straps and Velcro fasteners. First, draw the arm out of the sling and then from the knitted part of the support, then pull it downwards over the hand.

Material composition

Polyamide (PA), elastodiene/latex (LA), viscose (CV)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Close the Velcro fastening to avoid damage to other laundry items.
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

Medical device | Single patient – multiple use

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