# BORT OmoFX Vario

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**REF** 121 020

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D121020|2022-04|006 ML | Rev.01

### DE BORT OmoFX Vario

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### EN BORT OmoFX Vario

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined abduction position.

#### Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of impingement syndrome, rotator cuff tear, endoprosthesis implantation

#### Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material.

#### Fitting instructions

The BORT OmoFX Vario shoulder-arm abduction brace allows very comfortable fixation of the shoulder joint and elbow at an abduction angle of optionally 30° or 45°. Its removable cover ensures pleasant, hygienic wearing properties.

The OmoFX Vario can hold your arm in a 30° or 45° position. First, decide in which abduction position the arm is to be held. The arm is held in the 30° abduction position when the flatter inclination of the cushion faces upwards (towards the head). The arm is supported at the 45° abduction angle if the more inclined side of the cushion faces upwards (towards the head). The semi-circular side of the cushion is positioned to the side of the body. To effect a change of position, e.g. from 30° to 45°, remove the arm slings attached using Velcro (3 pieces), rotate the cushion 180° and re-attach 3 arm slings to the top of the cushion.

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Then place the side of the arm concerned on the foam cushion and use the three Velcro strips to position it firmly to the cushion. Following this, guide the strap attached at the upper arm level over the shoulders from the rear. Guide the strap back through the guide loop on the cushion and close it using the Velcro surface. Guide the straight waist strap around the back at waist height and through the guide loop and close it. If necessary, tighten both straps so that the brace is tight on the body.

#### Removal

To remove the OmoFX Vario shoulder brace, open all the Velcro fasteners, loosen all the Velcro fasteners in the reverse order, then remove the brace.

#### Material composition

PUR foam (PUR), polyamide (PA), cotton (CO), polyester (PES) The sewn in textile label provides the precise material composition

#### Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically Close the Velcro fastening to avoid damage to other laundry items. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 04.2022

Medical device | Single patient – multiple use

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