# BORT OmoStabil Desault-Weste

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
**REF** 121 200

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D12120012022-041004 ML 1 Rev.01

### DE BORT OmoStabil Desault-Weste

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### EN BORT OmoStabil Desault's Bandage

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined position.

#### Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of contusion, luxation, subcapital humerus fracture, rotator cuff lesion

#### Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.  
– remove the medical device prior to radiological examinations  
– if wearing at night is prescribed, avoid negative influence on the circulatory system  
– in case of numbness, loosen the medical device and remove it if necessary  
– in case of persistent complaints, consult the physician or specialist retailer  
– use the medical device in accordance with therapeutic needs  
– only use other products simultaneously after consultation with your physician  
– do not make any changes to the product  
– do not wear it on open wounds  
– do not use in case of intolerance of one of the materials used  
– no re-use – this medical aid is intended for treating one patient  
– whilst wearing the brace, please neither use any creams nor ointments on our around it as they can destroy the material

#### Fitting instructions

1 Depending on whether the vest should be worn on the left or right shoulder, turn it so that the semi-circular shaped shoulder pad (1) lies on the shoulder side to be immobilised and the lower arm support (2) hangs in front of the body.

2 As the support can be worn both on the right and the left, please cover the Velcro surfaces pointing outwards by turning the material. For easier handling, please attach strap 3 to strap 4 with Velcro. That simplifies the later attachment of strap 7 later.

3 Close the support in front of the chest by attaching the support part (6) to part (5) with Velcro.

4 Incline the forearm of the side affected (normally 90°), then close the lower arm support (2) towards the head.

5 Close the Velcro fastener for the elbow support (4) beneath the elbow so that it rests comfortably and securely in the designated recess.

6 Finally, close the remaining Velcro fastener (7), which is at upper arm level, at the front side of the body.

#### Removal

To remove the OmoStabil Desault Vest, open all the latching halves in the abdominal and arm areas, then remove the support.

#### Material composition

PES foam (PES), polyamide (PA)

The sewn in textile label provides the precise material composition.

#### Cleaning information

Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  
 Do not iron  Do not clean chemically  
Close the Velcro fastening to avoid damage to other laundry items.  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 04.2022

Medical device | Single patient – multiple use

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