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
**REF** 121 980

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D121980|2023-07|009 ML I Rev. 01

### DE BORT OmoSAT

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### EN BORT OmoSAT

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined position.

#### Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of distortion, contusion, subcapital humerus fracture, signs of paralysis

#### Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on our around it as they can destroy the material

#### Fitting instructions: The image depiction is from the right-hand side.

- Place the forearm in the arm pouch and position your elbow in the semi-circular stitched down area.
- Position the training ball in the metacarpus area and fix it using the Velcro area in the arm pouch.
- The arm pouch is closed in the direction of the antecubital fossa using the short forearm strap included. This strap fixes the forearm in the pouch.

4 Position the Velcro claw with the attachment module approx. 3 cm behind the wrist.

**Tip:** Here, the attachment module button should be in the direction of the body side.

5 The shoulder attached to the plastic ring with Velcro and the waist belt are guided jointly over the opposite shoulder and waist together in the threaded plastic buckle.

6 The plastic buckle is then hung in the interlock button. A light pull towards the hand region then engages the quick-release fastener and secures it.

**Tip:** The desired length setting and fixation can be set individually to the desired length, both using the shoulder and waist belt Velcro fasteners. This setting need only be made once.

7 The fitting process is completed.

8 **Tip:** You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

#### Removal

To remove the OmoSAT, move the plastic buckle towards the antecubital fossa and unhang the plastic buckle from the hanging button simultaneously. Open the Velcro fasteners in the reverse order and then remove the brace.

#### Material composition

Polyamide (PA), PUR foam (PUR), polyester (PES), cotton (CO) The sewn in textile label provides the precise material composition.

#### Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically Close the Velcro fastening to avoid damage to other laundry items. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 07.2023

**MD** Medical device |  Single patient – multiple use

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