BORT Immob-Schiene einteilig

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REF 145 070

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D14507012023-071005 ML | Rev.01

DE BORT Immob-Schiene einteilig

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EN BORT Immob Splint Single-Piece

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for the immobilisation of the knee joint in a defined (stretched/bent) position.

Indications

Pre-operative, post-operative, post-traumatic, muscle or ligament injuries.

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness or pressure sores, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

Fitting instructions

Before fitting the brace, the integrated aluminium splints located in the splint bags in the middle and side sections must be adapted to the patient's existing body contours.

Place the immobilisation splint around the leg from the rear so that the aluminium splints located at the back (dorsally) are in the middle of the popliteal fossa.

Guide the Velcro fasteners through plastic eyelets (D rings) and turn them to attach them with Velcro. Fit the splint so that despite immobilisation of the knee joint no congestion is caused in the calf and foot region.

Removal

To remove the device, loosen the belt straps and take the brace off the leg. Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition

Polyamide (PA), cotton (BW), polyester (PES)
The sewn in textile label provides the precise material composition.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items. Remove the aluminium splints before washing. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

Medical device | Single patient – multiple use

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