BORT Immob-Schiene einstellbar

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REF 145 080

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D145080I2022-10I003 ML I Rev. 01

DE BORT Immob-Schiene einstellbar

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a splint for the immobilisation of the knee joint in and adjustable position.

Indications

Pre-operative, post-operative, post-traumatic, muscle or ligament injuries.

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness or pressure sores, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

Information for technicians:

Important: The adjustable BORT immobilisation splint should always be adapted to the individual contour of the patient's leg by an expert cold forming (bending) the reinforcement splints.

Setting the desired flexion position: The splint is provided at an angle of 0° (full knee extension). If required, the joint angle can be changed to 15° flexion. For this, simultaneously remove the reinforcement splints on the rear of the legs from the splint bags and replace them with the light metal reinforcement splints provided in a bent position. Following this, remove the side reinforcement splints by loosening the Velcro connection and attach the bent reinforcement splints with a textile cover provided on both sides of the splint using Velcro.

Fitting instructions

Place the brace when sitting or lying from the rear around the stretched leg so that the kneecap is enclosed within the recess intended for this purpose.

Guide the Velcro straps through the guide loops and close the Velcro fasteners with medium tension so that a noticeable immobilisation of the knee joint is achieved, but no congestion is caused in the calf and foot area.

Removal

To remove the device, loosen the belt straps and take the brace off the leg.

Material composition

Polyamide (PA), cotton (CO), polyester (PES)
The sewn in textile label provides the precise material composition.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items and remove the aluminium splints. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2022

Medical device | Single patient – multiple use

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