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## BORT Stabilo Lumbar Back Brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### Intended purpose

This medical device is a brace for relief and functional mobilisation of the lumbar spine.

### Indications

Major lumboschialgia without absolute indication for surgery, when surgery is contraindicated, when surgery is rejected, major radicular and pseudo radicular lumbar syndrome (cautionary, therapy-resistant), spondylolisthesis with severe lumboschialgia, severe lumbar deformity with hypermobility associated with facet syndrome/osteoarthritis, condition after spine surgery, e.g. dorso-ventral/dorso-lateral fusion, spinal column compression (single level with deformity and multi-level), thoracic surgery with severe prolapse/lumbar instability, spinal canal stenosis with paresis, post-traumatic vertebral fractures with significant spinal damage.

### Contraindications

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

### Fitting instructions

#### Putting the brace on:

1. Open both the front clasps.
2. Unlock both side-positioned, round Boa fasteners by pulling them upwards.
3. Ensure that a gap of approx. 6 cm is given on the back between the pelvic girdle semi-shells. Should this not be the case, pull some of the pulley out of the Boa fasteners and pull the pelvic girdles apart.

4. Place the brace around the torso so that the upper parts of the plastic elements are positioned above the iliac crest.
5. Close the elastic support using the integrated hand straps.
6. Check the length of the grey straps, they may have to be tightened and shortened.
7. Close the lower grey fastener, then the upper grey fastener and check the correct fit of the brace. The upper part of the plastic shell should be positioned well above the iliac crest.
8. Then lock the Boa fasteners by pressing the rotary knob in. If you hear a "click", the Boa is locked. Stepless circumferential stabilisation can be set by turning it. You then feel the clear, circumferential support of the plastic stabilisation elements.

#### Adaptation to the course of therapy

The BORT Stabilo lumbar back brace can be adapted to the course of therapy, the decision regarding this is taken by the treating physician.

- Relief of the lumbar spine can be adjusted steplessly (see above).
- On the dorsal side, the lordosis reduction grade can be set by bending the aluminium splints which are integrated in the pelvic girdle.
- During the course of therapy, it may become necessary to reduce the stabilising effect of the lumbar brace. For this, remove the more rigid version of the splints from the externally attached splint bags on the lumbar side of the body support. The more flexible splints can then be inserted.
- If stabilisation is no longer indicated by the plastic elements, they can be removed from the carrier support which can then be worn as a normal support.

#### Removal

To remove the brace, unlock the Boa fasteners by pulling them out. To loosen the brace, the pulley from the Boa fastener can be pulled out using the black flaps. Open the grey straps and the front fastening and take the brace off.

#### Material composition

Basic support: Polyamide (PA), polyurethane (PU), polyester (PES), elastane (EL)

Pelvic girdle: Polyamide (PA), polyester (PES), polyurethane (PU)

The sewn in textile label provides the precise material composition.

#### Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Use a delicate wash for the grey support. Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items. Stretch back into shape and dry in the fresh air. Avoid heat.

Clean the grey pelvic girdles with a damp cloth.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before

submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

### **Useful life/Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 07.2023

**MD** Medical device |  Single patient – multiple use

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