BORT Stabilo FA Rückenorthese

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
REF 180 500

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D180500I2023-07I008 ML I Rev. 01

DE BORT Stabilo FA Rückenorthese

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EN BORT Stabilo FA Back Brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for partial relief due to lordosis reduction of the lumbar spine.

Indications

Facet joint arthrosis, intervertebral disc protrusions/prolapse, severe low back pain, lumboschialgia, lumbago, lumbago with sacroiliac joint involvement, chronic pain in the lumbosacral area, spondylolisthesis, muscular imbalance, segment instability, lumbar spinal canal stenosis.

Contraindications

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body, latex allergy.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

1 Before putting the BORT Stabilo FA Back brace on, open the front fastener. Loosen the 4 narrower tension straps so that it can be put on without restrictions. For this, undo the Velcro ends of the straps and place them loosely on the velour of the strap.

2 Place the BORT Stabilo FA Back brace around the waist and close the front fastener, using the ergonomic hand straps. Ensure that the back frame is placed centrally.

3 Tighten the two lower tension straps first.

4 Then tighten the two upper tension straps. Please regulate the tension of the tension straps in accordance with your needs and the shown indications.

Tip: The crocodile hand straps make it possible for you to remove them as required and to shorten the double velour individually to the requirements of the circumference measurements. Then, reattach the crocodile hand straps to the double velour.

Removal

To remove the brace, open the tension straps and the front fastening and take the brace off.

Material composition

Polyamide (PA), elastodiene/latex (LA), viscose (CV)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

Medical device | Single patient – multiple use

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