BORT StabiloPro Knieorthese

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
REF 182 150

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D18215012024-081008 ML I Rev.01

DE BORT StabiloPro Knieorthese

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EN BORT StabiloPro Knee Brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for supporting the physiological guidance of the knee joint with circumferential Velcro fastener restraints and side guidance splints with joints which are adjustable in their scope of movement.

Indications

Mild to moderate side ligament instability of the knee joint, meniscus injuries, gonarthrosis.

Contraindications

Thrombotic risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

The brace can be worn directly on the skin or over the clothes. Pay attention to avoid creasing the clothes in the vicinity of the brace as this can result in pressure sores.

Important note: The brace may only be adapted by expert staff.

Before putting the brace on, the side joint splints must be adapted to the existing contour of the leg by an expert. To adapt the aluminium splints, remove them from the splint bags. You can access the splint bags from the upper area of the brace, shown by the blue highlighted Velcro fastener. Shape the splints to the lateral contours of the leg using cold-shaping (if necessary with a bending iron). Finally, re-insert the joint splints into the joint pockets. Here, the blue extension restrictions must point towards the patella recess. **Tip:** For fitting the brace, first open the Velcro straps, guide the foot through the brace and then slip it over the leg affected. Position the brace so that the patella recess is under the kneecap. Then, close the Velcro fasteners under tension, avoid strangulation here.

Setting the extension/flexion limit: If required, it can be necessary to restrict the scope of the brace's movement. The extension is factory set to 0°, flexion is unrestricted.

1 Remove the outer joint cover. For this, lift the joint cover flap and slide the joint cover away downwards.

2 Important: The extension/flexion limits are inserted into the joint cover and not directly into the joint body.

The joints and the limit stop sets are marked with "A" and "B". Remove the existing limit stop from the joint cover and take the required limit stop from the limit stop sets provided, the angular degrees are noted on the limit stop wedges. Extension wedges (blue) are marked with an "E", flexion wedges (grey) with an "F". Insert the required limit stop into the joint cover and press it in place.

Attention: Be certain to store the inserts with the labelling 0° (blue) and 140° (grey) and insert them back into the joint cover after free movement has been achieved again. Both inserts serve to cover the round joint cover openings.

3 + 4 Finally, re-attach the joint cover. For this, attach the joint cover into the holder in the lower joint region and apply light pressure downwards until it engages. Depending on the limit stop you have selected, you must align the bending angle of the joint accordingly for joint cover assembly.

Example: You have inserted a 10° extension stop into the joint cover, therefore the joint must be bent slightly for assembly of the cover. The following restrictions are available:

– Extension: 0°, 5°, 10°, 15°, 20°, 25°, 30°
– Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Removal

To remove the brace, open the Velcro straps, grasp the lower edge of the brace at the side on the splints and pull it downwards over the foot.

Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.


Material composition

Polyamide (PA), polyurethane (PUR), elastane (EL), viscose (CV), cotton (CO)
The sewn in textile label provides the precise material composition.

 The product contains latex and can trigger allergic reactions.

Cleaning information

Before washing, remove the joint splints and close the Velcro fastener in order to avoid damaging other items of clothing.

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron  Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

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Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 08.2024

 Medical device |  Single patient – multiple use

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