# BORT CoxaPro Hüftorthese

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**Med**  
Mitten im  
gesunden  
Leben

## Sprachen Languages

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	03
		Anziehanleitung für den Anwender .....	04
		Anleitung für den Techniker.....	07
<b>EN</b>	English	Instructions for use .....	13
		Fitting instructions for the user.....	14
		Instructions for the technician.....	17
<b>FR</b>	Français	Mode d'emploi .....	23
		Instruction d'application destinée à l'utilisateur.....	24
		Instructions destinées au technicien.....	27
<b>ES</b>	Español	Instrucciones de uso .....	33
		Instrucciones de colocación para el usuario.....	34
		Instrucciones para el técnico.....	37
<b>IT</b>	Italiano	Istruzioni per l'uso .....	43
		Istruzioni per indossare l'articolo.....	44
		Istruzioni per il tecnico.....	47
<b>NL</b>	Nederlands	Gebraiksaanwijzing .....	53
		Aanwijzingen voor het aantrekken voor de gebruiker.....	54
		Aanwijzingen voor de orthopedisch instrumentmaker.....	57
<b>CS</b>	Čeština	Návod k použití .....	63
		Návod k navlečení pro uživatele.....	64
		Návod pro technika.....	67

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

**Intended purpose**

This device is a brace for the mobilisation of the hip joint with adjustable restricted movement in one plane for limiting and stabilising the movement of the hip joint when flexing and extending it.

**Indications**

All indications for which functional security and/or mobilisation of the hip joint is necessary, such as: condition after hip joint luxation, loosening of hip endoprostheses, medium to severe coxarthrosis, muscular imbalance of the lumbar pelvic hip region, e.g.: tractus iliotibialis syndrome, piriformis syndrome

**Contraindications**


Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

**Application risks/Important notes** ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

## Fitting instructions for the user

### Positioning the pads:


#### Trochanter pad

Place the trochanter pad provided on the intended Velcro surface in the inner area of the brace (see image).

#### Piriformis pad accessory REF 911 126 (not included in the scope of delivery)


If a trochanter pad is already attached in the intended area, you must first remove it by loosening the Velcro fastener cautiously.

Attach the piriformis pad to the area intended using Velcro (see image). Ensure that the pad is aligned correctly and the marked areas tally with the shape of the pad.


**Wearing information:** Both the trochanter and piriformis pads should exercise a noticeable but pleasant pressure on the hip region after the brace is fitted. The personal wearing period is determined in accordance with the instructions from the treating physician.

1


Place the opened brace around the hip and then centrally on the pelvis. The upper edge (sewn-in label) may neither be higher than the upper iliac crest nor extend to the waist. Push your middle and index fingers from the right and left hands into the flaps intended for this purpose on the ventral plate fastener. Then use both hands to pull the latch plates forwards together evenly. Place the left latching half on the stomach and Velcro the right latching half onto it. After attaching the fingers using Velcro, remove the fingers first from the left and then the right flap.

2


Ensure that the brace joint is positioned on the side of the thigh at the height of the hip joint ②, that is to say approx. 3–4 cm (approx. 2 finger's width) above the so-called trochanter major on the thigh bone. Due to its prominent shape, it can easily be felt ①.

3


Take hold of the non-elastic, narrow straps coming from the back and re-attach them with a pleasantly tight tension to the stomach fastener using Velcro.

4


Attach the upper narrow elastic restraint to the front area of the pelvic restraint with Velcro for easier fitting

**5**


Observe that, when closing both elastic restraints, the legs are positioned neutrally to the direction of movement during the further fitting process.

**6**


Take hold of the affected side of the hip using the designated finger flap on the side of the brace.

**7a**


**7b**


Then, use your other hand to guide your fingers into the guide flap of the wide lower restraint and pull it evenly around the thigh to avoid the formation of creases. Close the Velcro fastener flat on the outer side of the brace casing. If necessary, smoothen it again. The wide strap can also be temporarily attached to the knitted fabric Velcro fastener to simplify attachment.

**8**

Then, loosen the previously attached narrow elastic restraint from the pelvic strap and also guide it around the thigh crease-free under tension. Attach it again beneath the hip joint using Velcro. Close the Velcro fastener evenly. To simplify this, the finger flap can also be used here as required. Pay attention to tight fitting of the Velcro fasteners.

### Instructions for the removal of the pelvis buckle and joint before the washing process

Please take into account that the pelvis buckle ① with the joint ② must be removed prior to washing in order to avoid possible damages both to the items of clothing and the brace joint.

Start by removing the lower articulated shanks ② from the Velcro fastener area of the thigh cuff. Ensure at the same time that the position is marked in order to maintain the technician's original setting. In the next step, the narrow straps ③ which run in the pelvis buckle area are unthreaded from the guide loops ④.

After this, you can remove the guide flap ⑥ from the pelvis buckle ①. This guide flap is attached directly to the velour which is on the brace's knitted fabric.


For re-attachment of the pelvis buckle, carry out the steps in the opposite order.

The illustration shows the components already removed.


## Instructions for the technician

### Adaption of the course of the joint and splint

1

2


**1/2** Remove the pelvis buckle including the joint by detaching the lower articulated shank ② from the Velcro fastener area. Then, thread the narrow straps ③ which are guided around the pelvis buckle out of the guide loops ④ and passage loops ⑤ of the pelvis buckle.

Finally, you can remove the guide flap ⑥ from the pelvis buckle. Then, detach the pelvis buckle from the knitted fabric Velcro connection. Now, the individual components are freely accessible.


After the joint and pelvis buckle have been detached from the knitted fabric support, the side course of the splint can be adapted according to the anatomical conditions as follows.

3

Detach the screw-type connection which joins the brace joint to the pelvis buckle. This connection is located under the cover on the inside of the pelvis buckle in the region of the joint.


4

Then, move the joint horizontally on the pelvis buckle. This enables ideal adaptation of the brace joint to the anatomic joint. Once the desired position has been reached, re-tighten the screws.


**Information:** All Allen screws with an M4 thread should be tightened to a maximum torque of 1.5 Nm. Use the Allen key for screws with an across flat of 2.5 mm.

5

**5/6** When adapting the joint splints, pay attention that the course of the joint is in the patient's walking direction. Adaptation of the course of the joint both at sagittal and frontal levels may be necessary. This can be done using rounded bending hooks or manually with so-called cold forming.

6


## Setting extension / flexion without the use of tools


Check the markings on the joint cap and ensure that you use the respective markings for the right-hand or left-hand side of the hip joint.

For setting on the right, refer to the **R** marking and use the **L** marking for the left-hand side.

**1**


**2**


When adapting for the **RIGHT-HAND** side of the hip, the restriction for flexion can be found under the marking “FLE” on the right-hand side and the restriction for extension under the marking “EXT” on the left-hand side of the joint.

When adapting for the **LEFT-HAND** side of the hip, the restriction for flexion can be found under the marking “FLE” on the left-hand side and the restriction for extension under the marking “EXT” on the right-hand side of the joint.

**3a**


**3b**


For setting flexion and extension, use the wedge with the respective marking “A” or “B”. This marking is equivalent to the marking on the joint cap. Place it beneath the marking for flexion “FLE” in order to restrict flexion and under the marking “EXT” to restrict extension.

4


Flexion and extension positions can be restricted to 0°, 15°, 30°, 45°, 60°, 75° and 90°. The wedges are pushed into the lateral interlocking and locked downwards by being lodged. The A/B labelling of the wedges here is in your direction. Ensure that the angled wedge is set at the desired number of degrees.

5

Additionally, the joint can be freely positioned on the velour surface using Velcro. As a result, it is possible to react individually to the positioning of the hip joint (e.g. flexion contractures)

**Instructions for the removal of the pelvis buckle and joint:**

In as far as it becomes necessary within the scope of recommendations for the course of medical therapy to downgrade the brace whereby the pelvis buckle ① and the joint ② have to be removed from the textile part, you can do this as follows:

Start by removing the lower articulated shank ② from the Velcro fastener area. Then, thread the narrow straps ③ which are guided around the pelvis buckle out of the guide loops ④ and passage loops ⑤ of the pelvis buckle.

After this, you can remove the guide flap ⑥ from the pelvis buckle. Then, remove the pelvis buckle from the knitted fabric Velcro connection. The guide flap ⑥ is then attached directly to the velour which is on the brace's knitted fabric.

In order to use the brace without the pelvis buckle and joint, thread the narrow straps ③ through the guide loops ④ on the guide flap ⑥ towards the “bulging sail” on the front fastener. This ensures that the brace also sits correctly and comfortably when downgraded.

The illustration shows the components already removed.

## Removal

To remove the brace, please loosen the circular restraint around the thigh and open the frontal fasteners around the hip. Following this, you can remove the brace.

## Material composition


Cotton (CO), Polyamide (PA), elastodiene/latex (ED), elastane (EL), viscose (CV)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

## Cleaning information

 Delicate wash 30 °C  Do not bleach (only wash with colour or mild detergents)  Do not dry in a tumble dryer  Do not iron  Do not clean chemically

Before washing, remove the pelvis buckle with the joint splint and close the Velcro fasteners in order to avoid damaging other clothes items. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 01.2024

**MD** Medical device |  Single patient – multiple use

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