BORT Peroneus-Schiene

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REF 470 620

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D470620|2023-07|006 ML I Rev. 01

DE BORT Peroneus-Schiene

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EN BORT Peroneal Splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a peroneal spring, this brace holds the foot in its individual functional range in case of a tibialis anterior muscle paralysis.

Indications

Flaccid paralysis of the tibialis anterior muscle, reduced mobility of the patient, for example in the case of peroneal paralysis, lumbar nerve root damage, polio, stroke.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Protect the open injury on the foot with a bandage or similar. Always wear the brace together with a sock or tights and not on bare skin. The BORT peroneal splint can only be worn together with a shoe. Position the brace in a shoe and put it on to the foot affected together with the shoe. Pay attention to pressure-free fitting. Then, close the fastener on the lower leg.

Removal

First, open the fastener on the lower leg. Then, the fastener on your shoe.

Now, slip your foot out of the brace. Finally, you can then remove the brace from the shoe.

After removing the brace, check your foot and calf for erythema and pressure sores.

Tip: Close the Velcro fastener in order to avoid early wear to the Velcro surface.

Material composition

Polyurethane (PUR) 75 %, polyethylene (PE) 12 %, polyamide (PA) 13 %

Cleaning information

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. Close the Velcro fastenings to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

Medical device | Single patient – multiple use

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