IT

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EN

The treatment with Peroneal Splint serves to correct the leg positioning in case of the foot lifting and peroneal muscle paralysis, which untreated could result in pes equinus.

Instruction to use:

Protect open injuries on the foot with a dressing or similar. Do not wear the brace on the bare skin but always together with a stocking or sock. The BORT long soled Peroneal Splint can only be worn together with a shoe. Position the brace in the shoe and put the shoe together with the brace on the affected foot. Ensure a pressure-free fit. Now close the fastener on the lower leg. After removing the splint, check the foot and calf for skin reddening and pressure points.

Indications:

Flaccid paralysis of the tibialis anterior muscle, reduced mobility of the patient, for example in the case of peroneal paralysis, lumbar nerve root damage, polio, stroke.

Contraindications:

Peripheral arterial occlusive disease (PAOD), lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid, sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body.

Composition:

Plastic: 100 % polypropylene
 Knit loop/Velcro closure: 100 % polyamide
 Calf pad: 50 % polyamide, 50 % PU-foam

Important Note:

No reutilisation – this medical device is meant for the care of one patient only. Wear brace only when under strain, not during rest or sleep. Avoid hindering blood circulation. In case of numbness loosen brace and remove it if necessary. If complaints continue, please contact the doctor.

If pressure points, abrasions or similar occur, have your orthopaedic technician or doctor check the fit of the splint immediately.

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