BORT Knopfloch PIP Schiene

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REF 470 700

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D47070012022-041003 ML 1 Rev.01

DE BORT Knopfloch PIP Schiene

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EN BORT Buttonhole PIP splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for immobilisation of the PIP joint.

Indications

Buttonhole deformity (Boutonnière Deformity), immobilisation of the PIP joints.

Contraindications

Unclear soft tissue swelling, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Put the splint over the finger affected so that the PIP joint can be covered by the brace and the DIP joint can still be bent.

Please observe that the distal end of the brace must not reach beyond the distal antecubital fold.

To attach the splint, you can use a skin-friendly, self-adhesive fixing tape (not included in the scope of delivery) and place it circumferentially around the rear end of the splint, thus attaching the brace to the finger.

If the splint is too long, it can be processed using a funnel grinding machine at low rotational speed or shaped thermoplastically using heat (approx. 150°C). Only trained expert staff should carry out this processing. Do not shape the brace on the skin when heated. Avoid burns.

Removal

To remove the splint, remove the fastener used and pull the splint off cautiously.

Material composition

100% polyethylene (PE)

Cleaning information

- Hand-washing Do not bleach Do not dry in a tumble dryer
- Do not iron Do not clean chemically
- Do not use fabric conditioner. Dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/conformity

Status: 04.2022

Medical device | Single patient – multiple use

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