BORT Armtragegurt Mitella

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REF 610 120 CH

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D610120CHI2025-01I004 ML 1 Rev.01

DE BORT Armtragegurt Mitella

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EN BORT arm support strap Mitella

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is an arm support strap for immobilising the shoulder joint in a defined position.

Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of periarticular irritations and diseases in the glenohumeral joint region, simple fractures, supraspinatus syndrome.

Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid congestion of the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- when wearing the arm support strap: please neither use any creams nor ointments locally in the region of the attached aid as they can destroy the material

Fitting instructions

Guide the arm affected through the loop with the non-adjustable length until it is positioned beneath the elbow joint. Ensure that the strap is not twisted. The sewn arm strap position should be positioned towards the body. Then guide the strap behind the back and let it fall downwards vertically over the healthy shoulder. Following this, hang the arm affected into the adjustable loop from the wrist.

To change the arm position, you can open the silver fastener and adjust the length of the loop.



Removal

To remove the arm support strap, remove your hand from the adjustable loop. Then remove the strap from the lower arm cautiously and pull the arm support strap off.

Material composition

Polypropylene (PP): 100 %

Cleaning information

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron  Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 01.2025

 Medical device |  Single patient – multiple use

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