BORT Valco SOFT Hallux-Valgus-Schiene

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REF 930 030

4 005862 083471

D930030I2023-10|007 ML I Rev. 01

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DE

BORT Valco SOFT Hallux-Valgus-Schiene

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EN BORT Valco SOFT Hallux-Valgus-Splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a hallux valgus correction brace for correction of an x-malposition of the large toe. An adjustable fastener enables stepless, individual setting of the correction pressure.

Indications

Hallux valgus after hallux surgery to ensure the success of treatment, post-operative capsule relief.

Contraindications

Thrombosis risk, extreme varicosity, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material
- The brace can be worn both during the daytime and at night.
- The brace can be worn both on bandages and also in sufficiently spacious shoes.

Important:

The padded aluminium reinforcement splint should be shaped to the individual anatomical situation by bending before initial fitting in order to avoid pressure sores.

Fitting instructions

Open the brace casing fully in a first step and place it on the ground. Then, place the affected foot on the brace so that the side reinforcement splint ends in the front area at the tip of the large toe. With the next step, guide the sealing tape around the metatarsus and close the Velcro fastener using appropriate tension. Here, you must observe that the integrated splay foot pad should be placed as centrally as possible behind the metatarsus heads II - IV.

With the next step, guide the large toe tape around the large toe using appropriate tension and close the attached Velcro fastener.

If desired or prescribed by the doctor, toes II to V can be kept in an upright position using the enclosed toe correction loops. For this purpose, guide the loops around the respective toes and Velcro them on to the lower side of the support using appropriate tension. If dynamic care is required,

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95% Polyamid (PA), 3% Polyester (PES), 2% Elasthan (EL)

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remove the aluminium rod carefully from the side splint bag for the aid and replace it with the enclosed spiral spring rod. Ensure that the spiral spring rod is placed correctly. In order to avoid pressure sores in the region of the metatarsophalangeal joint of the toe, use the padding included in the scope of delivery. Place it carefully at the corresponding location on the metatarsophalangeal joint of the toe in order to ensure pleasant and comfortable use. The enclosed heel strap serves to attach the brace over the heel securely and is attached to the shell of the brace. It enables optimum securing of the position and guarantees that the brace sits firmly.

Removal

To remove the BORT Valco Hallux valgus splint, open the Velcro fastener of the brace casing and the large toe tape.

Tip: Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition

95% polyamide (PA), 3% polyester (PES), 2% elastane (EL)

Cleaning information

Waschen delicate wash 30 °C ☒ Do not bleach ☒ Do not dry in a tumble dryer ☒ Do not iron ☒ Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items. Please remove the aluminium splint and the splay foot pad from the respective pre-fabricated slit of the support sheath prior to washing it.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2023

MD Medical device | Single patient – multiple use

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