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**REF** 930 140

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D930140I2022-04I006 ML I Rev.01

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### EN BORT Forefoot Relief Shoe

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a therapeutic shoe/relief shoe for relieving the forefoot.

#### Indications

Post-operative, wounds/injuries of the forefoot area on the sole of the foot, post-operative following corrective osteotomy in the forefoot area.

#### Contraindications

No contraindications are known if used correctly.

#### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- you may not drive a vehicle whilst wearing the BORT forefoot relief shoe
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

##### Important:

Protect the open injury on the foot with a bandage or similar. Always wear a non-slip shoe with the same heel height on the healthy foot. Do not wear the shoe directly on the skin.

1. Please remove the pre-adapted toe cap completely. If required, it can be adapted as a cold or dust protector using a Velcro fastener on the shoe.
2. Open all the Velcro straps on the shoe fully, then open the shoe uppers' flaps.
3. Check the shoe for foreign bodies and remove them before putting it on.
4. Put your foot into the opened shoe.
5. Close the shoe uppers' flaps in order to avoid any strangulation. The flap above the heel can be re-adjusted.
6. Finally, close the Velcro tapes. They can be shortened, if required.
7. If required, you can use Velcro to finally adapt the toe cap to the forefoot area.

#### Removal


Open the Velcro fasteners in-full. If you have assembled the toe cap, fold it forwards. Then, open the shoe uppers' flaps. Remove your foot from the forefoot relief shoe.

**Tip:** Close the Velcro fastener in order to avoid early wear to the Velcro surface.

#### Material composition

Polyurethane (PUR), polyester, (PES), polyamide (PA)  
The sewn in textile label provides the precise material composition.

#### Cleaning information

 Hand-washing  Do not bleach  Do not dry in a tumble dryer  
 Do not iron  Do not clean chemically  
Do not use fabric conditioner. Hand wash the shoe using lukewarm water and mild detergents.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 04.2022

Medical device |  Single patient – multiple use

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