# BORT Verbandschuh Comfort

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**REF** 930 320/930 320 R/L

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D93032012024-101008 ML | Rev.01

### DE BORT Verbandschuh Comfort

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### EN BORT Comfort Orthopaedic Shoe

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a therapeutic shoe/orthopaedic shoe (short-term/long-term) and is mostly used for the treatment of acute damages to the foot and restores mobility in-part.

#### Indications

Post-operative and other wounds, even when treated with extensive and possibly padded bandaging, where the healing process takes place over an extended period of time and weight can be placed partially on the foot whilst walking, even outdoors. An additional insole is not required.

#### Contraindications

No contraindications are known if used correctly.

#### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- you may not drive a vehicle whilst wearing the BORT orthopaedic shoe
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

**Important:** For putting the orthopaedic shoes on and removing them, please use the zip on the heel. The Velcro fasteners on the orthopaedic shoe are only used to adjust the width and fit.

Protect the open injury on the foot with a bandage or similar. Always wear a non-slip shoe with the same heel height on the healthy foot. Do not wear the shoe directly on the skin.

Open the Velcro fasteners at the side and fold the flap forwards. As a second step, open the flap in the metatarsal. Then, the zip on the heel. Check the shoe for foreign bodies and remove them before putting it on. Place your foot

in the orthopaedic shoe. As the next step first close the zip on the heel, then the overlapping area over the metatarsal. Finally, position the central front flap centrally and fix it on both sides using Velcro fasteners.

#### Removal

Open the Velcro fasteners at the side and fold the flap forwards. As a second step, open the flap in the metatarsal. Finally, the zip on the heel. Remove your foot from the orthopaedic shoe.

**Tip:** Close the Velcro fastener in order to avoid early wear to the Velcro surface.

#### Material composition

50% Polyurethane (PUR), 35% polyester (PES), 15% polyamide (PA)

#### Cleaning information

Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron  Do not clean chemically  
Hand wash the shoe using lukewarm water and mild detergents.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

Medical device | Single patient – multiple use

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