# BORT Diabetiker-Verbandschuh

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
**REF** 930 400 R/L

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D93040012024-101007 ML 1 Rev.01

### DE BORT Diabetiker-Verbandschuh

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### EN BORT Diabetic Orthopaedic Shoe

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a therapeutic shoe/orthopaedic shoe (short-term/long-term) and is mostly used for the treatment of acute damages to the foot and restores mobility in-part.

#### Indications

Wounds in the foot and ankle region due to diabetes mellitus, peripheral artery occlusive disease (PAOD), further wounds in the foot and ankle region, post-traumatic, post-operative.

#### Contraindications

Extensive plantar wounds, no use without adequate bandaging/covering of the skin in the wound area.

#### Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- you may not drive a vehicle whilst wearing the BORT orthopaedic shoe
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

#### Fitting instructions

**Important:** The adaptation of the BORT Diabetic Orthopaedic Shoe must be made by an expert (orthopaedic technician, orthopaedic shoe technician)

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technician. Protect the open injury on the foot with a bandage or similar. Use crutches when the shoe. Always wear a non-slip shoe with the same heel height on the healthy foot.

#### Honeycomb insole in the BORT Diabetic Orthopaedic Shoe

Insole for plantar pressure relief. The individually removable honeycombs allow individual adjustment – post-operatively or for acute wound care. A thin surface layer fixes the processed area with Velcro fastening and hinders the adjacent honeycombs from slipping. Pressure is reduced.

- 1 Place a sterile film on the wound or the pressure ulcer and isolate the surface affected with a marker pen on the film.
- 2 The patient should then stand on the inner sole in the shoe, thus transferring the marking to the uppermost blue layer of the inner sole.
- 3 Carefully remove the honeycombs beneath the marking where the pressure relief should be provided. Attention, do not remove more than 7 honeycombs per pressure relief zone, otherwise targeted pressure relief and pressure redistribution is not guaranteed!
- 4 Glue the stabilisation plate (included in the scope of delivery) fully onto the sole with the honeycombs. The stabilisation plate is equipped with a self-adhesive film, thus fixing and stabilising the honeycombs around the relief zone and guaranteeing that it remains stable and is not changed.
- 5 Place the adapted honeycomb sole (with the stabilisation plate downwards) back into the shoe.

#### Attention:

We expressly recommend wearing bandage material or socks. If the product is not correctly handled, its function and optimum protection are no longer guaranteed. Extended hollows can result in fatigue symptoms and loss of function of the surface layer.

#### Removal

**Tipp:** Close the Velcro fastener in order to avoid early wear to the Velcro surface.

#### Material composition

50% Polyurethane (PUR), 35% polyester (PES), 15% polyamide (PA)

#### Cleaning information

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically

Hand wash the shoe and inner sole using lukewarm water and mild detergents. For this purpose, remove the inner sole and dry it at a warm, breezy location. Avoid direct heat exposure here.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical

device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

 Medical device |  Single patient – multiple use

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