

REF 137 120

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### EN Hammer Toe Correction Pad

**Intended purpose**  
This medical device is a device for the correction of dislocation of one or more hammer toes.

**Indications**  
Hammer toe dislocation

**Contraindications**  
Do not use in case of open wounds. In case of diabetes, inflammations or vascular diseases, consult a physician prior to use!

cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

**Useful life/Lifetime of the product**  
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

**Duty of notification**  
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for

**Application risks/Important notes**   
The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- > remove the medical device prior to radiological examinations
- > remove it during longer periods of rest, e.g. when sleeping
- > in case of numbness, loosen the medical device and remove it if necessary
- > in case of persistent complaints, consult the physician or specialist retailer
- > use the medical device in accordance with therapeutic needs
- > only use other products simultaneously after consultation with your physician
- > do not make any changes to the product
- > do not wear it on open wounds

your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact)

**Disposal**  
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

**Declaration of conformity**  
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2022

Medical device | Single patient - multiple use

- > do not use in case of intolerance of one of the materials used
- > no re-use – this medical aid is intended for treating one patient
- > whilst wearing the product, please neither use any creams nor ointments on or around it as they can destroy the material

**Fitting instructions**  
Fix affected toes on the side covered with textile using the fleece tie, tighten it on the reverse side and attach it to the Velcro. Protruding fleece tie can be shortened. For an optimum fit, wear socks or tights over the hammer toe correction pad if possible.

**Removal**  
Loosen the fleece tie on the lower side. Then, pull the pad away to the front over the toes!

**Material composition**  
Sealing tape: 91% polyamide (PA), 9% elastane (EL)  
Reinforcement: 100% polyethylene (PE)  
Micro Velcro: 100% polypropylene (PP)

**Cleaning information**  
 Hand-washing Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

**Guarantee**  
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or

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