
10 187 160 | 2022-10 | 0022 ML | Rev. 01

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EN Hallux valgus splint

Intended purpose
This medical device is a hallux valgus correction brace for correction of a valgus malposition of the large toe. An adjustable fastener enables stepless, individual setting of the correction pressure.

Indications
Hallux valgus malposition

Contraindications
Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes
The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use - this medical aid is intended for treating one patient
- when wearing the device: please neither use any creams nor ointments locally on or around it as they can destroy the material.

- the device can be worn both during the daytime and at night
- the device can be worn in footwear

- Fitting instructions**
- Open the PediSoft Textiline Hallux Valgus splint fasteners fully.
 - In the next step, pull the brace over the large toe (hallux). Here, the relief pad should be precisely on the same level as the metatarsophalangeal joint.
 - Close the Velcro fastener of the PediSoft Textiline Hallux Valgus splint in the metatarsal with an even pull. Pay attention to pleasant wearing comfort and avoid strangulation.
 - Pull the laterally positioned corrective support towards the heel with a well-dosed pull, here a visible correction of the hallux valgus (large toe position discrepancy) should take place.

The wearing time can vary from patient to patient and is determined by the existing anatomical circumstances. Remove the brace at regular intervals and check the skin appearance for redness or pressure sores.

Removal
To remove the PediSoft Textiline Hallux Valgus splint, open the Velcro fastener of the brace casing and the lateral corrective tape.
Tip: Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition
Polyamide (PA), polyester (PES), elastane (EL), cotton (BW)

Cleaning information
 Delicate wash 30°C Do not bleach Do not iron Do not clean chemically. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can

find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact

Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity
Status: 10.2022

Medical device | Single patient - multiple use

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