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
Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for relieving and guiding the knee joint via the 3-point principle with harnesses and adjustable joints for restricting the scope of movement.

Indications

Partial relief of the medial or lateral compartment in case of level II-III gonarthrosis, bone bruise (post-traumatic bone marrow edema), medial/lateral gonarthrosis as a result of a tibial plateau fracture.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

– whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Instructions for use for patients

The brace can be worn directly on the skin or over the clothes. Avoid easing the clothes in the vicinity of the brace, this may result in pressure sores.

- 1 Open all the brace's fasteners - harnesses as well as support sheaths.
- 2 Sit down and slip forwards to the edge of the chair. Bend your knee joint slightly (approx. 30°).
- 3 Place your leg in the brace. Ensure correct positioning. The joint should be positioned at the height of the centre of the patella (middle of the kneecap) and on the side of the leg.
- 4+5 Close the fasteners in the following order, always starting from the middle of the knee:

1. Support sheath fastener close to the knee on the lower leg
2. Support sheath fastener close to the knee on the thigh
3. Support sheath fastener away from the knee on the lower leg
4. Support sheath fastener away from the knee on the thigh
5. Circular strap away from the knee on the lower leg
6. Circular strap away from the knee on the thigh

Whilst tightening the harnesses, smoothen them in the tension direction with your free hand. This avoids twisting the brace on the leg.

Information: should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician.

Allow yourself time to become accustomed to the brace and increase the wearing time gradually. You can remove the brace during longer periods of rest. If you have questions regarding the use and wearing time, please refer to your physician or technician.

Removal

To remove the aid, open the harnesses and then take the brace off the leg. Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

Information for the technician:

Important note: The brace may only be adapted by expert staff. The BORT Generation Soft OA knee brace can provide for axial differences of up to 5°.

Application of model A:

- Use in case of varus dislocation/relief of the medial compartment on the right leg (the brace is positioned laterally).
- Use in case of varus dislocation/relief of the lateral compartment on the left leg (the brace is positioned medially).

Application of model B:

- Use in case of varus dislocation/relief of the lateral compartment on the right leg (the brace is positioned medially).
- Use in case of varus dislocation/relief of the medial compartment on the left leg (the brace is positioned laterally).

Adaption of the BORT Generation Soft OA knee brace on the patient:

- 1 Open all the brace's fasteners - harnesses as well as support sheaths.
- 2 The patient should sit down and slip forwards to the edge of the chair, the knee joint is bent slightly (approx. 30°).
- 3 Place your leg in the brace. Ensure correct positioning. The joint should be positioned at the height of the centre of the patella (middle of the kneecap) and on the side of the leg.

4+5 Close the fasteners in the following order, always starting from the middle of the knee:

1. Support sheath fastener close to the knee on the lower leg
2. Support sheath fastener close to the knee on the thigh
3. Support sheath fastener away from the knee on the lower leg
4. Support sheath fastener away from the knee on the thigh
5. Circumferential harness away from the knee on the lower leg
6. Circumferential harness away from the knee on the thigh

Whilst tightening the harnesses, smoothen them in the tension direction with your free hand. This avoids twisting the brace on the leg.

6 Use the enclosed key to set the required correction pressure when the brace is positioned, in case of a slightly bent knee joint and when the patient is sitting. For this purpose, adjustment screws are provided above and beneath the joint. Avoid overcorrection!

Important: Let the patient try walking. If required, the fitting correction pressure can be readjusted.

7 Setting the extension/flexion limit:

- If required, it can be necessary to restrict the scope of the brace's movement. The treating physician provides information regarding any necessary restrictions.
- Take the side joint splint out of the sheath fastener and remove the condyle pad from the joint.
- Remove the screw for the limit stop to be set.
- Remove the limit stop from the joint and insert the required angular degree.
- Fix the stop used with the screw provided, re-attach the condyle pad with Velcro and tighten the joint splint again in the sheath fastener.

The following restrictions are available:

- Extension: 0°, 5°, 10°, 15°, 20° and 25°
- Flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

Then, instruct the patient in the use of the aid.

Material composition

Sheath fastener: 59 % polyamide (PA), 10 % elastane (EL), 31 % PUR foam (PUR)

Fastening straps: 100 % polyamide (PA)

Padding: 60 % EVA foam (EVA), 40 % polyester (PES)

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link:
www.bort.com/conformity

Status: 05.2023

MD Medical device |  Single patient – multiple use

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